BureDouble 2.0

SystemRoMedic™

Instructions for use – English

Bruksanvisning – Svenska

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
Max Patient Weight: 150 kg / 330 lbs

Table of contents

English	3
Svenska	14
Norsk	25
Danske	36
Suomi	47
Deutsch	58
Netherlands	69
Français	80
EMC	90
Dimensions	94
Parts list	96
Product overview	100

The Bure Double 2.0 is equipped with twin electric drives for raising and lowering the Walker and widening the walking frame.

1. Intended use

Walkers are to be used as indoor, movable and portable assistive aids, which together with approved accessories are intended to be used for transfer from sitting to standing positions, and walking training. Walkers are to be used by persons with reduced muscular strength, who shall be able to stand and perform simple leg movements. Walkers are to be used for leaning on when walking and standing during early mobilisation training. Walkers and accessories are intended for indoor use only in Home Healthcare environment and Professional Healthcare facility environment.

2. Contraindications, Precautions, Warnings

Contraindications

- The device may not be used by patients above the maximum weight indicated on the label of the device.

Precautions

- Check that the device is used on a dry and flat surface.
- Check the device is correctly mounted/assembled before its first use.
- Check the device after every folding/disassembly after any transport.
- Check height adjustment and leg spreading motions and inspect the actuators or gas pistons full ranges.
- Activate the brakes whenever the device is not in use.
- Check the function of the brakes before use.
- Inspect the device regularly to detect any signs of damage. If there are any visible signs of damage, contact Direct Healthcare Group Customer Service.
- Check that the hand control does not show signs of wear before every use.
- Check that the hand control markings are in accordance with the lifting functions.
- Check the battery status of an electrically operated model before using the device.
- Warranty applies only if repairs or alterations are made by personnel who are authorized by Direct Healthcare Group.
- Ensure that there are no obstacles or people in the way of the device, when moving.
- Handle batteries with care. Do not drop.
- Use only batteries and cables that are intended for the device as instructed by the manufacturer.
- Verify that potential accessories that might be used are properly attached to the device before use.
- Low speed is recommended when moving with the device.
- Take care not to drive the device over thresholds with high speed or force. Approach the threshold with lower speed and communicate the coming threshold to the user.
- In case of an electrical model of the device, check that the hand control does not show signs of wear before every use
- Keep the instructions for use together with the product so that it can be used as a reference if necessary.

- Beware of moving parts during usage.
- Perform a risk assessment on the user before using the Walker to ensure they can stand and perform simple movements.

Warnings

- The device may not be used by infants, children or adults who are 146 cm or shorter.
- The caregiver must be able to read and understand the Manual/IFU of the product.
- Caregiver: Never use the device if you have not had a proper training on how to use the device including the accessories.
- Maintenance/Service staff: Never maintain the device if you have not had a proper training on how to use and maintain the device including the accessories.
- It is important to use only approved accessories to prevent unintended detachment of components and subsequently a fall that may lead to patient injury.
- Use careful and gentle maneuvers when moving the device.
- Perform maintenance/service of the device, according to the instructions in the IFU (Instruction for use), at least once every 12 months.
- No parts of the device shall be serviced, while the device is in use
- Accessories must be properly fitted and tested in relation to the user's needs and functional ability.
- Special care must be taken when using strong electrical power sources such as diathermy and the like so that diathermy cables are not placed on or near the device. In case of doubt consult with a Direct Healthcare Group representative.
- When lifting the device, only use the denoted handholds. At least two people are required for lifting the device.
- Do not leave a patient unattended when the device is being used.
- The device can not be used by patient alone.
- Never move the walker by pulling on the actuator device or gas piston.
- The device must not be lowered into water.
- The device must not be cleaned using steam.
- The device must not be left or stored in a damp or humid environment.
- The device must not be charged in a wet room.
- The device must not be used outdoors, only indoors.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective grounding.
- No modification of this equipment is allowed
- Do not modify this equipment without authorisation of the manufacturer.
- If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of this equipment.
- Do not use the device while it is charging.
- Do not use the device if any of the cables shows signs of wear and tear.
- Do not put the device in direct sunlight as it may get hot.
- Do not run over the cables with the walker or other devices.

Always read the user manual

Always make sure that you have the correct version of the manual. The most recent version of all manuals are available for downloading at/from our website: www.directhealthcaregroup.com.

Before first use:

- Inspect the packaging for any damage.
- Check for correct product being delivered.
- Check to ensure that the device is free from damage.
- Check if all the parts/components are included in the package.
- Check if all approved accessories are included in the package.
- Check if Product Quality Approval document is included in the packaging. Save this document for future contact with manufacturer.
- Make certain the device is properly assembled.
- Check lifting function and bottom-frame adjustment function.
- Check if it is possible to move the device in all directions.


3. Using Bure Walker

Standing up and walk training

1. Ask the user to sit on the edge of the bed or chair. Bring out the walking table and place it in front of the user. Lock the brakes on the wheels.
2. Ask the user to place their arms on the armrests. If necessary, ask the user to move forward on the seat base so that the user reaches the handles.
3. Ask the user to lean forward. To start raising the walking table, press the UP button on the controller. Raise the walking table to a desired height.
4. Unlock the brakes. Adjust the height of the walking table to a desired position
5. Assist the user while walk training.

4. Using Bure Walker with Gait Training Kit 56-388

- The support harness has four blue straps with a plastic buckle allowing for adjustments in length, having steel carabiner clips at the ends to be fasten to the gait walker. The carabiner clips on to the metal mounts placed under the table of the Gait walker.
- Apply the walker support harness by clipping on the two front carabiners to the front mounters, one on each side of the table.
- Place the user standing up in the walker. Pull back the harness in between the legs and clip on the carabiners to the back mounters, one on each side of the table. Adjust the harness to the user by pulling in the back strap loops.
- **Observe!** Intended use of the trainer support harness is to practise walking. It is **not** a sitting or lifting aid.

Standing Plate

Attach standing plate by hooks on the main frame of the walker.

Raising/Lowering the Frame

The desired height is adjusted with Arrow <UP> or Arrow <DOWN>. Adjust bottom frame using <WIDER> arrow or <NARROWER> arrow

Handset Plus has charge indicator / battery (low, medium, high) and Service indicator (see further under "Hand Control Plus").

Adjusting the Handles

To adjust handles to the desired position; loosen screw, adjust handle, tighten screw.

Left and Right Adjustment of Armrests

The armrests can be adjusted sideways for optimal support; loosen screw, adjust armrest, tighten screw.

NOTE: The rear attachments on the frame/armrest are interchangeable (left with right and vice versa) to provide different widths. The locating pins on the cushion's underside/rear can also be moved. These adjustment capabilities make customisation possible.

Using the Parking Brakes on Castor Wheels

When parking the device, activate the parking brake on the lockable castor wheels. Lock the wheels by pressing the green pedal brake down with your foot. To release the brake, press the green protruding part with your foot. The caregiver should assist the patient in activating the parking brake in case the patient cannot activate them him or herself.

Using Hand Brakes

Running brakes: Pull the brake levers towards you to reduce speed.

You will brake the rear wheel that is located on the same side as the brake lever.

If you apply only one of the brakes, the walker will turn to that side. When the lever is released, the braking will stop. Do not use the brakes continuously while the walker is in motion.

Parking brakes: Push the brake levers forwards from you to activate the parking brakes.

Be aware that the lever has some resistance before it stops and locks the rear wheel.

Pull the brake levers back towards you to release the parking brakes.


Hand Control Plus

**1. Service indicator – flashes red after 11 months' operation.

Constant light (red) after 12 months' operation.

The indicator is a reminder to review the maintenance instruction (see Maintenance Instructions). This should be done at least once every 12 months.

Reset: hold down buttons <2a> + <2b> for at least 30 seconds. **The indicator will go out around 30 seconds after the hand control is activated (battery conservation function).

Battery Charging

Battery charging must take place when the device is not in use.

Charging starts automatically when power plug is connected to 230V wall socket. Green light on the control box will turn on, once the power plug is connected.

- The LED flashes green during charging.
- ON/OFF: To use the product, press ON.
- The green lamp will begin to flash.
- Use the OFF button (hold down for 5s) to switch off the battery if the product will not be used for an extended period.
- The LED will show a constant green light when the battery is fully charged.
- When the battery is in normal state/in operation, the LED flashes green at REGULAR intervals, namely 0.5s ON, 4s OFF.
- When the battery requires charging (remaining charge <25%) the LED will flash orange/red at REGULAR intervals, namely 0.5s ON, 4s OFF.

NOTE: The battery has an in-built device that automatically switches it off when remaining charge is <20%. The LED is then extinguished. The battery must now be fully charged before the electronics can be used again.

Control box

The control box cover can be opened by using a screwdriver or similar; refer to the instructions which are sitting on the control box.

NOTE: Make sure that plugs are connected to the correct sockets otherwise the actuators can be damaged.

1. Socket for the height adjustment actuator.
2. Socket for the width adjustment actuator.
3. Hand control socket.
4. Battery terminal socket.
5. Plug socket (NOTE: Plugs must ALWAYS be connected to the control box for IP classification to be valid).

Trouble shooting

If the lifting function cannot be activated, check the following:

- That all cables are properly and securely connected, pull out the contact and plug it in again firmly.
- That battery charging is not in progress.
- That the battery is charged.

If the device is not working properly, contact your distributor.

If the device makes unusual noises:

- Try to determine the source of the sound. Take the device out of operation and contact your distributor.

Accessories

Art. No.	DESCRIPTION
58-328-1	Heel strap Standing plate
56-336-4	Oxygen holder
56-338	Drip rod complete with attachment
56-337	Side support, adjustable
56-337-PU	Side support, adjustable PUR
56-352	Lower Leg Support
56-388	Gait training belt
57-368	Guide handle
56-389	Soft basket with attachment
56-333-DB	Standing plate Short Bure DB, Rise DB
56-333/SW-DB	Standing plate Short, SW, Bure DB, Rise DB
56-383-DB	Standing plate Long Bure DB, Rise DB
56-383/SW-DB	Standing plate Long, SW, Bure DB, Rise DB
56-339/75-2-0	Manual brake, 75mm wheel, 2 handles
56-339/75-2-1	Manual brake, 75mm wheel, 1 handle
56-339-2-0	Manual brake, 2 handles
56-339-2-1	Manual brake, 1 handle
56-305	Touch-up paint Bure, Silver 8
56-306	Touch-up paint Bure, Dark grey RAL 7021
58-375-X2	Hand control PLUS TIMotion TH10

5. Maintenance

The device must undergo thorough inspection at least once per year. Inspection must be performed by authorised personnel and in accordance with Direct Healthcare Group's instructions.

Repairs and maintenance may only be done by authorised personnel using original spare parts.


Check that	Action if a problem is noticed
<p>General The Walker feels firm/free from play. The Walker does not rattle during manoeuvres. The Walker is level and all of its castors are in contact with the floor. The Walker is not dirty.</p>	<p>Clean with lukewarm soap solution or alcohol-based cleaning agent (no petroleum products).</p>
<p>Armrests Armrests are intact and clean. Armrest widening is functional.</p>	<p>Clean; replace if damaged. Replace lock knob. Fit lock screws/plates as required.</p>
<p>Handle Grips are not damaged/dirty. Handle adjustment is functional.</p>	<p>Clean with lukewarm soap solution (no petroleum products). Replace handle grips. Fit lock screws/plates as required.</p>
<p>Electronics Check that the charging cable is connected to the control box. Check that the hand control is connected to the control box. Check that the battery, actuators and control box are not loose.</p>	<p>If any part needs to be replaced refer to the wiring diagram in the control box wiring section. Fit new charging cable; this must always be connected to the control box. Connect or fit new hand control as required. Tighten; replace with new fasteners as required.</p>
<p>Height adjustment Raising and lowering are functional. The Walker feels stable at maximum settings. There is no play between the vertical frame and bottom frame. Height adjustment clamps lock. The electric actuator fitting must be free from play. Height adjustment via hand control is functional. The electric actuator must run smoothly at a constant speed.</p>	<p>Fit new guide sleeves (in frame) or end plugs (in chrome tubes) as required. Tighten all lock bolts. Replace clamps. Tighten the attachment fitting concerned; replace bolts/lock nuts as necessary. Check that the battery/hand control/ actuator are connected in accordance as in the control box wiring section. Charge the battery. The part list shows the design and which fasteners must be checked.</p>
<p>Frame parts There is no mechanical damage. There are no scratches. End plugs/lower frame fitted</p>	<p>If mechanical damage is present, contact DHG Customer Services. Touch up as necessary. Fit new end plugs.</p>
<p>Castors/brakes Castors roll easily/tread not damaged. The castors are firmly fastened to the lower frame. Castor brakes functional on all castors. Handbrake functional.</p>	<p>Clean or replace castors. NOTE! Castors are always fixed to the lower frame with thread locker or lock nuts. The castors are sealed and we do not recommend their dismantling; instead, replace the whole castor assembly. Tighten the castor bolt and fit locking nut or use thread locker (depends on model). Replace castor assembly. Adjust the brake or fit new castors.</p>

Power cord may only be replaced by DHG service personnel or technical personnel trained by DHG.

DHG can provide upon request necessary technical information in order to assist maintenance and repair of the device.


Storage and transportation

The walker should be transported and stored in 5°C to 50°C and in normal humidity, 20% to 80% non-condensing. The air pressure should be between 700 and 1060 hPa. Leftmost symbol indicates storage and transportation. Let the walker reach room temperature before the batteries are charged or the walker is used. The walker should not be stored so that it is exposed to dust, or so that the battery or gas piston is exposed to direct sunlight.

Operation

The operating environment should be 5°C to 40°C, relative humidity 20% to 80% non-condensing, and atmospheric pressure 700 to 1060 hPa. See also Technical specification below. Leftmost symbol indicates operating condition.

Cleaning and disinfection

After every use:

Clean the walker armrests and handles. Wipe using a mild soap solution/surface disinfectant. However, the product is not corrosion proof.

Weekly and possibly more often depending on frequency of use:

If necessary, clean the whole walker with a cloth with warm water or a soap solution and check that the castors are free from dirt and hair. The electronic components (battery, control box, hand control, actuators, cables) should be cleaned with a damp cloth only. To avoid degreasing of the piston rods, the actuators should be retracted to minimum stroke and without load before cleaning. Ensure that the device is dried thoroughly after cleaning.

If disinfection is needed:

70% ethanol, 45% isopropanol or similar should be used.

Do not steam clean due to risk of corrosion.

Do not use cleaning agents containing phenol or chlorine, as this could damage the materials.

With environments with high risk of infections use the walker with PU-cushions (Accessory number 57-361-L and 57-361-R)

General Maintenance for gait training Sling:

Read the product label.

Do not use rinsing agent.

To prolong product life, avoid tumble-drying.

Recycling/Disposal

When the product is no longer usable it should be recycled in accordance with legislation and regulations in the country concerned. All electrical parts including batteries must be removed and recycled as electrical components. Ask your local recycling station for further information about how the different materials in the product (metals, plastics, electronics) should be recycled. Used batteries can also be returned to Direct Healthcare Group or a Direct Healthcare Group distributor for recycling.

6. Technical Specifications

Control box	TiMotion Voltage in 220V-240VAC, Voltage out 24V, 50Hz, 1250mA, 70VA, Encapsulation IP66, Electrical safety Class II (Double insulated) Duty cycle: 10%, Max. 2min. ON/18 min./OFF
Batteries	TiMotion Rechargeable Li-Ion Battery, 24 V - 4A, 96Wh, Encapsulation IP66, Capacity 4,0Ah, Cable length 600mm. Duty cycle: 10%, Max. 2min. ON/18 min./OFF
Lifting Actuator	TiMotion 24VDC, Max 8,0A, Push: 2000N, Pull: 2000N, Self-locking capacity (min): 1000N, Current: 6.6 (+/-1,0)A, Stroke 400mm +/-2mm, Mounting dim. 605mm, Encapsulation IP66. Duty cycle: 10%, Max. 2min. ON/18min. OFF
Widening Actuator	TiMotion 24 VDC, Max 8,0A, Push: 1000N, Pull: 1000N Self-locking capacity (min): 1000N, Current: 2.6 (+/-0,4)A, Stroke 55mm +/-2 mm, Mounting dim. 167mm, Encapsulation UIP. Duty cycle: 10%, Max. 2min.ON/18min. OFF
Hand Control	Cable length straight 300mm + coiled 400 mm, IP66
IP class	IPX5
Material	Steel
Castors	125/100/75mm
Weight	33kg
Expected Service life	8 years
Max patient weight	150 kg / 330 lbs
Operating environment	5 °C to 40°C, relative humidity 20% to 80% non-condensing, and atmospheric pressure 700 to 1060 hPa.
Storage and transportation environment	5 °C to 50 °C, relative humidity 20% to 80% non-condensing, and atmospheric. Environment pressure 700 to 1060 hPa.
Materials in contact with human skin	Materials that user and operator/caregiver are exposed to during usage of the product for prolonged time: Poroflex (synthetic leather cushions) P ^x Polyurethane (on PU cushions) P Dryflex (handles) P PA6 + 20% GF (brake lever) P Hand control C ^x
Applied parts:	Armrest P Handles P Knee pad P Foot rest support P Foot rest support clamps C Handles adjusting knobs C Hand controller and its holder. C Handles to adjust the armrest width. C Hand control C Charger C Charger cable C Footplate

^x P - patient, ^x C - caregiver

7. Product label


Bure DB Timotion, 125mm MD CE

SystemRoMedic

REF 56-303

(01)07331769043117(11)220301(21)610231377


U In: 220-240 VAC, 50 Hz
I In: Max 1250 mA





IPX5
Direct Healthcare Group Sverige AB
Torshamnsgatan 35, 16440 Kista, SWEDEN
www.directhealthcaregroup.com
Made in Sweden

UDI (Unique Device Identification) can be found within the barcode human readable text. It contains.
(01) EAN code
(11) Date of Manufacture. The date is in the format YYMMDD
Y = Year, M = Month and D = Date.
(21) Serial Number displayed as five digit code.

Rev.02

8. Table of Symbols

	This product complies with the requirements of the Medical Device Regulation 2017/745		Visual Inspection
	Medical Device		Read the manual
	Caution		Product Code
	Legal manufacturer Direct Healthcare Group Sverige AB		Batch Code
	Type BF applied parts, according to the degree of protection against electric shock		The device is intended for indoor use
	WEEE Symbol May not be discarded in domestic waste		Class II Equipment
	Duty Cycle: 2 min in active (ON) mode. 18 min in rest (OFF) mode.		Maximum patient weight
	Walker Height		Walker Width
	Do not step on the device		Refer to instruction manual (IFU)
		Weight (mass) of the device, the Maximum Patient Weight, and total sum. All in kg.	

9. How to report a serious incident

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the MHRA or another competent authority of the country in which the user and/or patient is established.

UK

info@directhealthcaregroup.com

T: +44 (0) 800 043 0881

F: +44 (0) 845 459 9832

Other:

info.export@directhealthcaregroup.com

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Electromagnetic Compatibility (EMC)

Electromagnetic emissions and test levels

The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Emission test	Standards	Compliance
RF emissions	CISPR 11	Group 1
RF emissions	CISPR 11	Class B
Harmonic current emissions	IEC 61000-3-2	Class A
Voltage fluctuations and Flicker emissions	IEC 61000-3-3	Complies

Electromagnetic Immunity

Acceptance criteria for the EMC pass criteria is unintentional movement above $\pm 10\text{mm}$ is not allowed. The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.

Acceptance criteria for the EMC pass criteria is unintentional movement above $\pm 10\text{mm}$ is not allowed. The product is intended to be used in the stated environments with electromagnetic levels as specified below. The caregiver and/or user of the product assure that the product is used in such an environment.


Radiated Fields in close proximity, Immunity test and compliance. Standard: IEC 61000-4-39			
Dwell time: 3 sec			
Test levels (A/m)	Modulation Pulse modulation	Mod. Frequency (kHz)	Test Frequency
8	CW	CW	30 kHz
65	50% Duty cycle	2.1	134.2 kHz
7.5	50% Duty cycle	50	13.56 MHz

Immunity test	Basic EMC standard or test method	IEC 60601-1-2 Edition 4 Test levels and compliance	
		Professional healthcare facility environment†	Home healthcare environment
Surge	IEC 61000-4-5	± 0.5 kV, ± 1 kV	
Voltage dips	IEC 61000-4-11	0% U _r ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
		0% U _r ; 1 cycle and 70% U _r ; 25/30 cycles. Single phase: at 0°	
Voltage Interruptions	IEC 61000-4-11	0% U _r ; 250 / 300 cycle	
Conducted disturbances induced by RF fields	IEC 61000-4-6	6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	6V in ISM and amateur bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
		Dwell time: 3 sec Frequency step size: 1%	
Radiated RF Electro-magnetic fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz
		Dwell time: 3 sec Frequency step size: 1%	
Electrical fast transients / burst	IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	
RATED Power frequency magnetic field	IEC 61000-4-8	30 A/m 50 Hz or 60Hz	
Electrostatic Discharge (ESD)	IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV Number of discharges: 10 for each polarity	

Radiated Radio – Frequency, Proximity fields from wireless communication equipment & compliance. (Dwell time: 3 sec.)
 Standard: IEC 61000-4-3


	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test levels (V/m)
Spot Frequencies	385	380 to 390	TETRA 400	Pulse modulation 18 Hz	27
	450	430 to 470	FMRS 460, FRS 460	FM ± 5 kHz Deviation 1 kHz sine	28
	710	704 to 787	LTE Band 13, 17	Pulse modulation 217 Hz	9
	745				
	780				
	810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
	870				
	930				
	1720	1700 to 1990	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 217 Hz	28
	1845				
	1970				
	2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
	5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
	5500				
5785					
Other Identified frequencies	433	-	-	Pulse modulation 2 Hz	3
	-	860-960	-	Pulse modulation 2 Hz	54
	2450	-	-	Pulse modulation 2 Hz	54

Dimensions


	A	B	B	C	D
Wheel diameter	Turning diameter	Armrest height	Armrest height PU	Wheel area diameter	Length
125	1214	900-1300	870-1270	218	950
100	1186	880-1280	850-1250	190	935
75	1184	850-1250	820-1220	188	930

Parts


No.	Art No.	Description	No.	Art No.	Description
1	58-366	Middle frame	28	57-787	Washer 10x25x1,5
2	-	Widening assembly Bure Double, R&G, see separate view	29	57-885	Screw MF6S M6x10
3	58-335	Telescopic Inner pipe	30	57-811	Serrated lock washer M16
4	58-365	Top frame	31	57-803	Locking washer M10 6.1x11.8x1.6
5	-	Included in 58-366	32	57-816	Screw MC6S M10x40
6	56-340 57-310	Guide sleeve Middle frame bushing	33	57-804	Washer 22x6,5x8
7	58-379-L	Bottom beam left Bure double	34	57-808	Pin bolt M6x25
8	58-379-R	Bottom beam Right Bure double	35	57-792	Screw MC6S M8x35
9	57-306	PU Armrest cushion Right	36	57-806	Washer 8.4x16x1.5
10	57-305	Armrest cushion Right	37	57-812	Tapping screw RXS 2,9x4,5
11	57-049 57-319 57-051	Castor wheel 125mm brake Castor wheel 100mm brake Castor wheel 75mm brake	38	57-769	Washer 12,5x35x2
12	57-775	Dome cover wheel	39	57-766	Lock nut CL6 M12
13	57-776	Dome cover wheel top	40	57-879	Cap nut M6
14	58-360	Cover Wheel attachment	41	57-864	Screw K6S M6x20
15	57-041-3	Holder remote control Bure double	42	57-884	Screw MC6S M10x45
16	58-373/ Lion	Control box 220V Digital Li-On TiMotion Bure	43	58-354-1	Wheel attachment
17	58-389	Battery box 4.0 Ah TBB5/Li-On, see separate view	44	57-801	Lock nut M8
18	56-353	Bure Handle assy	45	58-326	Spiral Cord EUR - T-Motion, Bure
19	57-919	Star knob handle	46	58-374-X2	Hand control TiMotion TH7
20	57-917	Wing nut	47	58-377	Actuator Timotion 2000N/5200rpm Bure
21	58-338	Bushing middle frame	48	57-904	Washer 6.4x12.5x2 nylon
22	56-308	Cable holder Bure charger	49	57-828	Pop rivet AL/ST 4,8x10
23	57-639	Screw K6S M6x8	50	57-369	Battery clamp Alu
24	57-805	Washer BRB 6,4x12x1,6	51	57-370	Battery bracket Alu
25	57-817	Lock nut CL6 M10	52	57-883	Screw MF6S M6x40
26	57-825	Washer 10,5x18x2	53	57-878	Screw K6S M6x25
27	57-848	Nordlock 10,7x21x2,5			

No.	Art No.	Description
17		Battery box 4.0 Ah TBB5/Li-On

No.	Art No.	Description
2		Widening assembly Bure Double, R&G

No.	Art No.	Description	No.	Art No.	Description
1	58-356	Housing Turn motor	14	57-826	Washer 8,4x29x1,5
2	58-352	Crossection bottom Bure Double	15	57-760	PVC Spacer 8x12x6,4
3	58-378	Actuator TA16 TiMotion	16	57-765	Screw MF6S M8x25
4	58-364	Turning wheel	17	58-302	Crossection bottom bushing
5	57-763	Ball joint angled M8	18	57-818	Screw K6S M6x10
6	58-381-1/2	Threaded steel rod M8	19	57-778	Washer Nylon 05-8-12-10
7	58-380-R	Turning tube right	20	57-774	Washer Nylon 13x34x1,5
8	58-380-L	Turning tube left	21	57-876	Nut M6M M8
9	57-801	Lock nut M8	22	57-882	Screw K6S 10.9 M8x55
10	57-806	Washer 8.4x16x1.5	23	57-767	Washer Starloc 8x15,3
11	58-358	Cover Housing Turn motor	24	57-771	Screw RXS 3,5x6,5
12	57-777	Washer nylon 05-8-12-3	25	57-875	Set screw SK6S M8x6
13	57-836	Spring lock washer M8			

Product Overview

REF

Article number	Product name
56-303	Bure DB Timotion, 125mm
56-303-AUS	Bure DB Timotion, 125mm, AUS
56-303-B/0-PU	Bure DB Timotion, 125 mm, handbrakes, PU
56-303-B/1-PU	Bure DB Timotion, 125 mm, 1 handbrake, PU
56-303-PLUS	Bure DB Timotion, 125mm, Plus
56-303-PLUS-PU	Bure DB Timotion, 125mm, PU, Plus
56-303-PU	Bure DB Timotion, 125mm, PU
56-303-PU-AUS	Bure DB Timotion, 125mm, PU, AUS
56-303-UK	Bure DB Timotion, 125mm, UK
56-303/100	Bure DB Timotion, 100mm
56-303/100-PU	Bure DB Timotion, 100mm, PU
56-303/75	Bure DB Timotion, 75mm
56-303/75-B/0	Bure DB Timotion, 75mm, handbrakes
56-303/75-B/0-PU	Bure DB Timotion, 75mm, handbrakes, PU
56-303/75-UK	Bure DB Timotion, 75 mm, UK
56-303/75-PU	Bure DB Timotion, 75mm, PU
56-303/75-PU-UK	Bure DB Timotion, 75mm, PU, UK

Main configuration is following: 125 mm wheels, no handbrakes, cushions made of synthetic leather, regular hand control, no special features. Following variations are possible:

- 125/100/75 mm – indicates diameter of wheels
- Handbrakes – indicates that the model has a handbrake on both handles
- 1 handbrake – indicates that the model has only 1 handbrake on one of the handles
- UK, AUS, US – models designed specifically for these countries
- Timotion- means that the model includes electronics from Timotion for adjustment of the height
- Plus – means that the model includes a Plus hand control
- PU – means that the model includes cushions made of Polyurethane
- RS – means that the model has rotation stop.

Simple solutions for great results

SystemRoMedic® is the name of Direct Healthcare Group's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic® is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic® is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic® offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website;
www.directhealthcaregroup.com

For questions about the product and its use

Please contact your local Direct Healthcare Group and SystemRoMedic® representative. A complete list of all our partners with their contact details can be found on our website; www.directhealthcaregroup.com.

Moving Health Forward

CE Type BF Applied Parts

Medical Device Class I. The product complies with the requirements of the Medical Devices Regulation 2017/745.

Direct Healthcare Group Sverige AB
Torshamnsgatan 35,
SE-164 40 Kista, Sweden
Tel: +46 (0)8-557 62 200
info@directhealthcaregroup.com
www.directhealthcaregroup.com