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Maintenance instructions for bathroom adaptations

Maintenance, washing and cleaning instructions

Washing, cleaning and maintaining aids is a necessity not only for hygiene and safety, but also for prevention of infections. Below, you will find general information on maintaining our products. For more information, read the instructions in the manuals. Our products are constantly being further developed and improved; therefore we reserve the right to modify products without prior notification. The most recent information is always available on our website: www.directhealthcaregroup.com/nl.

Maintenance before and during use

Direct Healthcare Group advises to inspect the product regularly on the following points:

- Make sure that the product is properly installed and/or assembled.
- Make sure all screws and bolts are tightened properly.
- Make sure clips are in the right position and can safely perform their duties.
- Regularly check that all the components are in place and properly tightened.
- Check for possible damage.

Daily cleaning

1. For optimal hygiene it is recommended to rinse the product immediately after use with hot water or a mild soap solution (PH 6-8). We recommend not to use chlorine, solvents, abrasive or aggressive cleaning agents, as this can cause discoloration of the coating, PUR and plastic parts, or even may cause damage the material itself. Some parts, such as a seat or bedpan can be removed or folded to facilitate easy cleaning. It is important that these parts are put back correctly.
2. After rinsing, the product should be wiped with a clean, damp, soft cloth. The rinse and wipe ensures that soap and secretions are removed immediately and so avoid the risk of bacteria growth.
3. Allow the product to dry, preferably in a well ventilated area, before it is used again.

Disinfection

The product should be disinfected according to the protocol of the facility with regard to method, duration and frequency. Below are the general guidelines from Direct Healthcare Group:

The products, including the polyurethane parts, can tolerate temperatures up to 85°C for three minutes. Alternate options are: 90°C for one minute or 80°C for ten minutes. The disinfectants that have been approved by the local healthcare authorities may be used for disinfection. A 70% ethanol solution, 45% isopropanol solution or something comparable can be used without problems. Reuse if the product is reused or used for a new client, defective or worn parts must be replaced. Disinfection of the product as described above is necessary.

 85°C 185°F	 3 min
 90°C 194°F	 1 min
 80°C 176°F	 10 min

Recycling

When the product has reached its technical lifespan, the product must be disposed of in parts. The plastic parts can be sorted as residual waste. Metal parts should be treated as metal waste for optimal recycling.

**Direct Healthcare
Group**

Advancing Movement & Health[®]