Carina350EM/Carina350EE

SystemRoMedic™

User manual - English

Bruksanvisning - Svenska

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

REF

60600009

60600011

60600012


60600014

SWL: 165 kg/350 lbs

Table of contents

| | |
|----------------------------|---------|
| English | 4 |
| Svenska | 15 |
| Norsk | 26 |
| Dansk | 37 |
| Suomi | 48 |
| Deutsch | 59 |
| Nederlands | 70 |
| Français | 81 |
| Italiano | 92 |
| Español | 103 |
| | |
| Dimensions | 114-116 |
| Serial number | 117 |

1. Lift arm
2. Mast
3. Handle bar
4. Battery
5. Emergency stop button on control box
6. Rear castors with brakes
7. Sling bar with safety latches
8. Front castors
9. Emergency lowering (manual)
10. Motor/actuator for lift arm
11. Locking handle
12. Pedal for manual width adjustment (EM)
13. Tension strap
14. Quick-connecting locking pin
15. Hand control

Carina350 is a foldable, mobile patient lift adapted for safe and easy lifting of users weighing up to 165 kg/350 lbs. Carina350 is available in four variants with manual and electrical base widening, respectively and with standard height or low legs for use with lower beds. All variants are sturdy, stable and safe and, at the same time, light-weight for a steel lift. Quite simply, Carina350 offers a unique combination of functionality, performance and quality in relation to price. Carina350 is supplied complete and ready to use fitted with a stable sling bar with reliable safety latches and a hand control.

Safety information

Visual inspection

Inspect lift functions regularly. Check to ensure that material is free from damage.

Before use

Make certain the lift is properly assembled.

Check sling bar connection and safety latch function.

Check lifting function and base-width adjustment.

Check to ensure that the quick-connecting locking pin for the actuator is correctly installed.

Always read the user manual

Always read the user manuals for all assistive devices used during a transfer.

Keep the user manual where it is accessible to users of the product.

Always make sure that you have the right version of the user manual.

The most recent editions of user manuals are available for downloading from our website, www.directhealthcaregroup.com.


Under no circumstances may the lift be used by persons who have not received instruction in the operation of the lift.

It is strongly prohibited to modify the original product.

Assembly

Check to ensure that the following components are included

- Mast with lift arm, lift motor, sling bar and control box.
- Undercarriage with base-width adjustment pedals.
- Quick-connecting locking pin for mounting actuator.
- Hand control and cord.
- Locking handle for mounting mast on undercarriage.
- User manual.
- Charging cord.


The lift is folded when delivered. Set up the lift, according to the following instructions (no tools required)

1. Lock castor brakes.
2. Loosen the locking handle at the base of the mast (image 1).
3. Pull the mast forward, lift it up and place it in the base.


The mast **MUST** be fully seated in the base. Mast and base shall be aligned (image 2). Improper fit could cause the parts to fail or break.

4. Tighten locking handle (image 1).
5. Loosen the tension strap that holds the lift arm and lift up the lift arm.
6. Remove the actuator from the mount on the mast and move it to the mount on the lift arm.

It is important the lift arm is fitted correctly in the actuator (image 3). Improper fit could cause the parts to fail or break and serious risk of injury.

7. Check all mounting fixtures and locking handles (images 1).
8. Deactivate the emergency stop (image 4) and perform a final inspection (see Final inspection).



To fold the lift, follow the instructions below

1. Lock castor brakes.
2. Run the lift arm to the lowest position, and then press in the emergency stop.
3. Remove the actuator from the lift arm.
4. Lower the lift arm and secure it to the link on the sling bar with the tension strap.


Warning for risk of pinching (between boom and mast).

5. Loosen the locking handle at the base of the mast (image 1). The locking handle doesn't need to be removed.
6. Lift up on the mast and lower it towards the undercarriage.
7. Pull the mast back and tighten the locking handle (image 1)


Final inspection

Inspect the lift for signs of wear and damage. Check all castors and castor brakes.

Check all connections and fixtures including screws and bolts. Check the packaging to ensure that there are no loose parts.

Check the emergency stop function by activating the emergency stop button, and then pressing either the up or down button. If nothing happens when the up or down buttons are pressed, the emergency stop is functioning properly.


Grasp the hand control, press the up button and run the lift arm all the way up. Then, press the down button and run the lift arm all the way down.

Test base-width adjustment function by pushing down on the respective pedals for widening and narrowing the base. For EE use the buttons on the hand control.

Test lift function by lifting a person (not a user) using an approved sling. At the same time, check the emergency lowering function with someone in the lift, see section on Emergency lowering.

If the lift is functioning correctly, connect the charger and check to ensure that the charging lamp on the control box lights up.

NOTE!

Before the lift is used for the first time:

- it must be charged for at least 4 hours. See section on charging batteries.
- the control box service counter needs to be reset. To reset the service counters press both lift buttons on the hand control at the same time for 5 seconds. An audio signal will indicate that the timer has been reset.

Keep the user manual where it is accessible to users of the product.

Using the product


Important Information

- The lift must be assembled according to the instructions provided by Direct Healthcare Group.
- The lift may only be used on a level floor.
- Lifting accessories must be properly fitted and tested in relation to the user's needs and functional ability.
- Do not leave the user unattended during a transfer situation. 
- Under no circumstances may the safe working load be exceeded. See section on Safe working load.
- Never move the lift by pulling on the actuator! 
Do not push
- The lift must not be used or immersed in or under water including showering.
- The lift should not be left / stored in a high humidity environment.
- The lift must not be charged in a wet room.
- To ensure optimal function, the lift must be inspected regularly. See section on other maintenance details.
- Warranty applies only if repairs or alterations are made to the product by personnel who are authorized by Direct Healthcare Group.
- Using lifting accessories other than those approved can entail a risk. 
- Note that there is a pinch point between the manual emergency lowering mechanism and the lift arm when the lift arm is in its highest and lowest positions.
- Verify that all lifting accessories hang vertically and can move freely.

Safe working load

Different products on the assembled lift system; ie. lift, sling bar, lifting sling and other accessories may have different safe working loads. It is always the lowest maximum load allowed indicated for the respective products on the system that applies for the entire system. Always check the markings on all products in the unit and contact Direct Healthcare Group with any questions.

Charging batteries

A tone sounding when using the lift indicates that the battery need recharging.
Charge the lift after use to ensure that the battery is always fully charged.
Lock the castor wheels when charging the battery.

1. Connect the charging cable to a power outlet and in the connection for charger cable.
The symbol for charging is shown on the display.
2. Check to ensure that the lamps on the control box light up. The green LED lamp indicates that the charger is receiving power and the yellow LED lamp indicates that the battery is charging.
3. Charging stops automatically when the battery is fully charged.

Wall-mounted charger


1. Remove the battery pack from the lift and place it in the wall-mounted charger.
2. Check to ensure that the LED lamp on the front of the charger lights up.

NOTE!

Before the lift is used for the first time, it must be charged for at least 4 hours.
For maximum battery life, charge batteries regularly. We recommend daily charging when the lift is used daily.
The emergency stop must be deactivated during charging.

Battery information on the display

The battery discharging will be shown in four stages:

Battery state 1: The battery is ok, no need for charging (100 - 50 %).

Battery state 2: Battery needs charging. (50 - 25 %)

Battery state 3: Battery needs charging. (Less than 25 %) A tone sounds when a button is pressed in this battery state.

Battery state 4: The battery needs charging. (17V or lower) At this stage some of the functionality of the lift is lost. At this battery stage it is only possible to drive the lift arm down. Furthermore an audio signal will sound when a control button is activated. The symbol will switch between the two pictures for 10 seconds.

The battery symbol is shown when the control box is active until power down (2 minutes after use). It is not possible to use other battery types than BAJ1/BAJ2.

The battery level is measured via voltage. This means that it is possible to experience e.g. that the battery switches from state 1 to state 2 and back to state 1.

| | |
|---------|---------|
| ↑ | 12034 |
| A+S | 1257000 |
| 1kg/10s | 7 |
| SELECO | 90/360 |

Service information read-out

Basic service information can be read out on the display. To get the service information on the display please press the lifting arm up button for half a second. The information will be shown for ½ minute or until other buttons are activated.


- Total cycles done by the actuator
- Total work done by the actuator (ampere usage times seconds in use)
- Total number of overloads
- Days since last service/Days between services

Hand control

Raising/lowering the lift arm

Symbol indicate direction of travel.

Motion stops as soon as the button is released.



Electrical base-width adjustment (EE)

Markings on the buttons of the hand control indicate function.


Motion stops as soon as the buttons are released.

Manual base-width adjustment (EM)

Push down on the respective pedals at the back of the lift to widen and narrow the base.

Changing the battery

Emergency stop



To activate the emergency stop

Press the red emergency stop button on the control box.

Resetting

Turn the button in the direction of the arrows until the button pops out.

To prevent battery discharge, we recommend that the emergency-stop button is pressed in when the lift is not in use.


Emergency lowering

Manual emergency lowering

For manual emergency lowering, turn the round plastic knob on the actuator clockwise.

Electrical emergency lowering

For electrical emergency lowering, use the down button on the control box. Use a narrow object such as a pen.

Trouble shooting

If the lift adjustment cannot be activated, check the following

- That the emergency stop button is not pressed in.
- That all cables are properly and securely connected. Pull out the contact and plug it in again firmly.
- That battery charging is not in progress.
- That the battery is charged.

If the lift is not working properly, contact your dealer.

If the lift makes unusual noises

- Try to determine the source of the sound. Take the lift out of operation and contact your dealer.

Accessories

End cover Carina

Covers the end of the mast when collapsed.

End cover, article no.: 70200030

ReadyStrap

Locking strap, article no.: 16010

Fixed actuator bolt

If the lift is not intended to be folded for storage the spring pin can be replaced with

Fixed actuator bolt, article no.: 70200091

Lifting slings

Direct Healthcare Group's SystemRoMedic™ line includes a wide range of functional and comfortable,

high-quality lifting slings that are adapted for all types of lifting and for users with different needs.

The lifting slings are available in several materials and in sizes ranging from XXS to XXL. There are also special lifting sling models in XXXL and XXXXL for extremely large and heavy users. All models are safe and very easy to use.

Sling bars

SlingBar is a two-point aluminium sling bar available in three variants with different width. All variants of SlingBar have safety latches which prevent the sling straps from creeping out of the sling bar and all variants are for users weighing up to 300 kg/660 lbs. To get more room in a lifting sling when using SlingBar two-point sling bar, SlingBarSpreader M side bars can be used as an accessory. SlingBarSpreader M open up the lifting sling and provide for a more reclined position.

SlingBar XS article no.: 70200071

SlingBar S article no.: 70200001

SlingBar M article no.: 70200002

SlingBar L article no.: 70200003

SlingBarSpreader M article no.: 70200042

Hand control

Hand control HB33-6 for 60600009 and 60600014 with service, battery status and overload indicators

Hand control HB33-6, article no.: 70200089

Maintenance

The lift must undergo thorough inspection at least once per year. Inspection must be performed by authorized personnel and in accordance with Direct Healthcare Group's service manual.

Repairs and maintenance may only be done by authorized personnel using original spare parts.


Used batteries are to be left at the nearest recycling station. Used batteries can also be returned to Direct Healthcare Group or a Direct Healthcare Group dealer for recycling.

Cleaning/disinfection

If necessary, clean the lift with warm water or a soap solution and check that the castors are free from dirt and hair. Do not use cleaning agents containing phenol or chlorine, as this could damage the aluminium and the plastic materials. If disinfection is needed, 70 % ethanol, 45 % isopropanol or similar should be used.

Storage and transportation

If the lift is not to be used for some time or e.g., during transport, we recommend that the emergency stop button be pressed in. The lift should be transported and stored in -10°C to $+50^{\circ}\text{C}$ and in normal humidity, 30% -75 %. The air pressure should be between 700 and 1060 hPa. Let the lift reach room temperature before the batteries are charged or the lift is used.

Service agreements

Direct Healthcare Group offers the possibility of service agreements for maintenance and regular testing of your mobile lift. Contact your local Direct Healthcare Group representative.

Symbols

| | |
|--|---|
|  | White/blue Read user manual |
|  | May not be discarded in domestic waste |
|  | The product complies with the requirements of the Medical Devices Directive 93/42/EEC |
|  | Type B, according to the degree of protection against electric shock |
|  | The device is intended for indoor use |
|  | Class II equipment |
|  | Important Be observant |
|  Do not push | Red/ black Do not push or pull the lift by the actuator |
|  | Yellow/Black Pinch warning |

EMC

This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the equipment/system or shielding the location.

Technical Information

| | |
|--|--|
| Lifting speed | 23 mm/s without load. |
| Batteries | Two 12 V, 2.9 Ah valve-regulated, sealed, lead accumulator (gel-type batteries) |
| Charger | Max. 400 mA |
| Motor (mast) | DC 24 V, 5 A. IP X4. Operationtime: 10% at maximum continous running of 2 minutes, maximum 5 switching cycles per minute. Push: 6000N. |
| Motor (base) | DC 24 V, 5 A, IP X4. Operationtime: 10% at maximum continous operation of 2 minutes, maximum 5 cycles per minute. Push: 3000N. |
| Sound level | With load: upwards: 43 dB(A) downwards: 44 dB(A). |
| Material | Powdercoated Steel |
| Emergency lowering | Manual and electrical |
| Castors | Front 4", 100 mm, back 4", 100mm |
| Weight | 33 kg/73 lbs |
| IP class | IP X4 |
| Expected lifetime | 10 years |
| Operating forces buttons on handset | 4 N |
| Max. load | 165 kg/ 350 lbs |


Medical Device Class I. The product complies with the requirements of the Medical Devices Directive 93/42/EEC.

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

Dimensions Carina350EE/EM (60600009/60600011)


All measurement are in mm/inch

Dimensions Carina350EM/EE, with low legs (60600012/60600014)



All measurement are in mm/inch

Dimensions

- Turning radius is 615 mm/24.2 inch
- The total weight of lift is 33 kg /73 lbs f
- The weight of the heaviest component is 17,5 kg / 38.58 lbs
- D is movement in forward direction

All measurements are in mm/inch

Serial number

SystemRoMedic™

Simple solutions for great results

SystemRoMedic™ is the name of Direct Healthcare Group's unique easy transfer concept, the market's widest and most comprehensive range of clever, easy-to-use and safe transfer and lifting aids designed to make life easier, both for the user and for the caregiver. SystemRoMedic™ is a complete solution that provides for the majority of patient transfer or manual handling requirements. From the simplest to the most complex scenarios, from the lightest to the heaviest. The concept encompasses assistive devices for four different categories of transfers:

- Transfer, assistive devices for manual transfers of users between two locations.
- Positioning, assistive devices for manual repositioning of users within the same location.
- Support, assistive devices for mobility support e.g., during sit-to-stand or gait training.
- Lifting, assistive devices for manual and mechanical lifting of users.

Improved work environment, improved quality of care and cost savings

The philosophy behind SystemRoMedic™ is focused on the prevention and reduction of occupational injuries while allowing users to experience a greater sense of independence and dignity. Through a unique combination of training and a complete range of efficient transfer aids, SystemRoMedic™ offers improvement of both work environment and quality of care and, at the same time, achieves significant cost savings.

Always make sure that you have the correct version of the manual

The most recent version of all manuals are available for downloading at/from our website;
www.directhealthcaregroup.com

For questions about the product and its use

Please contact your local Direct Healthcare Group and SystemRoMedic™ representative. A complete list of all our partners with their contact details can be found on our website; www.directhealthcaregroup.com.

Direct Healthcare Group

Advancing Movement & Health

Direct Healthcare Group Sverige AB
Torshamnsgatan 35,
SE-164 40 Kista, Sweden

Tel: +46 (0)8-557 62 200
info@directhealthcaregroup.com
www.directhealthcaregroup.com