# PositioningSling, disposable

Instructions for use - English

Bruksanvisning - Svenska

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


**REF**

547810

**SWL: 454 kg/1000lbs**

## Table of contents

<b>English</b> .....	<b>3</b>
<b>Svenska</b> .....	<b>8</b>
<b>Norsk</b> .....	<b>13</b>
<b>Dansk</b> .....	<b>18</b>
<b>Suomi</b> .....	<b>23</b>
<b>Deutsch</b> .....	<b>28</b>
<b>Nederlands</b> .....	<b>33</b>
<b>Français</b> .....	<b>38</b>
<b>Italiano</b> .....	<b>43</b>
<b>Español</b> .....	<b>48</b>
<b>Dimensions</b> .....	<b>53</b>

## Intended use

The PositioningSling, disposable is designed to allow a client to be lifted and repositioned in a supine manner. It may also be utilized to roll or turn a client who requires frequent change of position while in bed. The sling can also provide support for a person while (s)he is lying on his/her side. The PositioningSling, disposable can be used in combination with ceiling lifts or mobile lifts, and can assist in client transfer from bed, shower trolley/stretchers and even off the floor. It has twelve attachment points, six on both sides, which can be connected to the sling bar hooks.

The sling is made of disposable material for **single client use**, which is advantageous for Infection Control and Prevention in care environments.

## Indication of use

Thorough risk assessment of the client, the environment and the task to be performed is always required when deciding which sling type and application is to be used. Always refer to your policies and procedures to ensure appropriate techniques are being used.

PositioningSling, disposable has 12 lift straps, 6 pieces along both sides.

To simplify handling for the caregiver, the lift straps have three loop alternatives which are marked with colors: white, grey and red. This makes it easier to use the same strap loop at both sides of the client to lift level.


The PositioningSling, disposable is marked with a product label, which is positioned at the head of the bed with the label facing the surface. An arrow on the label marks the middle of the sling for placing it under the spine of the client.

## Check safety

### Visual inspection

Check the condition and function of the assistive device regularly. PositioningSling disposable is a client specific product and should be disposed of when it becomes dirty or soiled. It must not be used if it has been laundered as the product integrity will be compromised. Before use, check for cleanliness, check that the seams and material surface are free from damage and are fully intact. If any part of the PositioningSling disposable are worn, faded or damaged it must not be used, it should be discarded and a replacement PositioningSling, disposable used instead. Apply load to the device and check to ensure that the loops withstand heavy loads. When starting to lift, gradually increase the load, checking the security of the lifting loops and integrity of the attachment points and the position of the end user throughout.

## Always read before use

Always read the instructions for use for all assistive devices used, including the lift system, prior to starting a transfer. Keep the instructions of use where it is accessible to trained users of the products.

Always make sure that you have the right version of the instructions for use.

The most recent editions of instructions for use are available for downloading from our website, [www.handicare.com](http://www.handicare.com).


- Always read the instructions for use for all assistive devices used, including the lift system, BEFORE USE.
- Keep the instructions for use where it is accessible to users of the product.
- Do not leave the client unattended during a transfer.
- To prevent discomfort and the risk of the client sliding out of the sling, trial fit the sling carefully.
- The lowest allowable safe working load always determines the safe working load of the assembled system.
- Always check the safe working loads for the lift and accessories before use. Do not exceed safe weight capacity.
- Please consult with a health care professional trained in safe patient handling for more details or uses of this or other sling models. Only people who have received proper training should attempt to use our lift systems.
- Improper use can lead to serious injury.
- Do not use the lift systems for anything but their stated purposes.

## Use of Positioning Sling, disposable

---

### Applying the sling

1. Roll the client onto his/her side. Maximize her/his participation.
2. Fold the sling in half lengthwise with the label showing on the outside of the fold. Lay the folded edge along the length of the client approximately one hand width away from the client's spine/body.
3. Roll the top folded layer of the sling towards the client. Make sure all straps are tucked into the roll. Push the material down into the mattress and under the client.
4. Return the client to supine and roll her/him to the other side.
5. Unroll the top layer of the sling out until flat. Return the client to supine. The sling should be flat under the client with straps free on all sides.


### Sling bar orientation

The sling bar can be oriented widthwise or lengthwise to the client. However, widthwise orientation is recommended to ensure a balanced weight distribution during lifting. For turning of the body, for example for turning in bed, the sling bar can be oriented lengthwise.

### Attaching the straps to the sling bar - general applying information


All straps, positioned from the head to the knees, must be attached correctly (top 4 or 5 upper most straps) on both sides of the PositioningSling, disposable to the sling bar. The sling bar in combination with the lift needs to be approved by Handicare for use with the PositioningSling, disposable. It can be a 2-, 4- or 8 point bar.

The bottom straps (usually below the knees) can be attached at the client's discretion. When leaving the bottom straps disconnected, the caregiver is required to guide the lower extremities. Take care to eliminate as many friction points as possible

#### Lifting – attaching the sling symmetrical to the sling bar

Attach the straps to the sling bar. To maintain a leveled lifting position, use the same color loops on the left and right sides of the sling. Make sure to always check that the loops remain attached to the sling bar during lift or transfer. When the client is not lying in the middle of the sling, the caregiver may switch to sling loops of different lengths, as long as there is no tension on the loops. The caregiver should make sure to stretch the sling loops simultaneously and to check that the client is positioned horizontally.



#### Higher up in bed

Start raising the lift and move the client higher up in bed.

### Turning in bed, attaching the sling asymmetrical to the sling bar

If the client is to be positioned in the center of the bed after turning, it is recommended to first reposition the client (see: Higher up in bed) laterally to one side of the bed and then roll the client towards the center of the bed (1).

To perform the turn, connect the four center-most lift straps on one side of PositioningSling, disposable (side opposite the direction of the intended roll) to the sling bar (2). If bed rails are not used, a caregiver should be on the side you are rolling the client toward to ensure the safety of the client. Raise the lift to gently turn the client towards the center of the bed.

### Change of PositioningSling, disposable

When replacing an old PositioningSling, disposable with a new one, start with a lift as previously described. When the client is lifted clear of the bed, place the new PositioningSling, disposable centered on the bed.

Lower the client back down to the bed and disconnect the used PositioningSling, disposable from the sling bar. Connect the new PositioningSling, disposable to the sling bar and roll the client onto their side (see: Turning in bed) and roll the used PositioningSling, disposable behind the client, then roll the client back down onto the bed. Do another turn with the new PositioningSling, disposable in the opposite direction and pull away the used PositioningSling, disposable from behind the client, then roll the client back down onto the bed.

**Leave in place**

When repositioning is completed, tuck the straps in under the mattress.

Best practices are to have the PositioningSling, disposable left directly under the client at all times. A flat sheet may be placed between the client and the sling, but results of an independent lab study demonstrate that the PositioningSling is safe to lay on, even with specialty mattresses and fragile client skin (per facility approval). The breathability, temperature and humidity against client skin are comparable to a standard hospital sheet. Assess client regularly for skin health.

**Tip!**

---

- Use with plain sheets on top of PositioningSling, disposable. For easier positioning a slide sheet, for example WendyLett, can be recommended.
- Place the client in the center of PositioningSling, disposable.
- A balanced fit is a safe fit: position the sling under and around the individual.
- The extra material in length available with PositioningSling, disposable allows the client to be repositioned several times before you need to reposition the client on PositioningSling, disposable.
- Recommended for use with traversing (x-y axis) ceiling lift systems.
- If possible, consider moving the bed in the first place instead of draw the client when he / she is lifted from the surface. When using a straight rail overhead lift system, please consider the position of the bed underneath the ceiling rail so as to ensure the client does not move unexpectedly when lifted.
- When positioning is completed, the straps/loops of the PositioningSling, disposable can be tucked under the mattress for safety and neatness.
- For washing/dressing/medical care of the lower part of the client, a positioning pillow/wedge can be put behind the back of the client before lowering the sling bar. This makes it possible for the client to keep the side-lying position, without the sling being in the way for (self-)care activities.
- Best practices are to have the PositioningSling, disposable left directly under the client at all times.
- For comfort and hygiene purposes place a flat sheet between the client and the sling.
- Assess client regularly for skin health.

## Symbols

	Medical Device Class I. The product complies with the requirements of the Medical Devices Directive 93/42/EEC.
	Blue/white Read user manual
	Middle of the sling
	Visual inspection
	Mechanical inspection
	Do not wash
	Do not iron
	Do not bleach
	Do not tumble dry
	Do not dryclean

### Material

Polyester

### Care of the product

Read the product label.

Do not launder. Discard the product if it is soiled or if there is a need to use another product. Discard the sling if any tears or serious defect are found. Single client use only.

 Medical Device Class I

The product meets the requirements of the Medical Devices Directive 93/42/EEC.

Contact your local distributor if you have any questions about the product and its use. See [www.handicare.com](http://www.handicare.com), [www.prismmedical.ca](http://www.prismmedical.ca) and [www.prismmedicalinc.com](http://www.prismmedicalinc.com) for a complete list of distributors.

### Expected life time

This client specific sling should last 1 to 3 months, if not soiled or laundered.

All our slings are compatible with our complete range of sling bars. To ensure anyone using our lift systems is properly trained, an on-going training program should be established.

Contact your representative should you have further questions.

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
## Drawing

Sling dimensions	cm	inch
Width	132	52
Length	198	78
6 straps on both sides, length	40.6	16

## **Always make sure that you have the correct version of the user manual**

The most recent version of all user manuals are available for downloading at/from our website; [www.handicare.com](http://www.handicare.com).

## **For questions about the product and its use**

Please contact your local Handicare representative. A complete list of all our partners with their contact details can be found on our website; [www.handicare.com](http://www.handicare.com).

***Handicare offers solutions and support to increase the independence of disabled or elderly people as well as to improve the convenience of those who are caring for them.***

*The Handicare Group is one of the leading healthcare companies in Europe with its own manufacturing organizations and sales companies in Norway, Sweden, Denmark, Germany, the Netherlands, Great Britain, France, China, Canada and the USA. Handicare's products are also distributed by partners in more than 40 countries worldwide. Our wide range of high-quality products includes a complete easy transfer system and other patient handling aids, stairlifts, car adaptations and bathing and toileting products.*

**PRISM MEDICAL**  
A Handicare Company

Prism Medical Ltd.  
485 Millway Avenue Unit 2,  
Concord ON L4K 3V4 Canada  
1.877.304.5438 [www.prismmedical.ca](http://www.prismmedical.ca)

10888 Metro Court, Maryland Heights, MO  
63043 USA Toll Free 1.866.891.6502  
[www.prismmedicalinc.com](http://www.prismmedicalinc.com)

Handicare AB  
Torshamnsgatan 35  
SE-164 40 Kista, SWEDEN  
Tel: +46 (0)8-557 62 200  
Fax: +46 (0)8-557 62 299  
E-mail: [info@handicare.se](mailto:info@handicare.se)  
[www.handicare.com](http://www.handicare.com)