Shower stool round

Linido

Manual - English

Manual - Svenska

Manual - Danske

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Art.no.: 200160001

SWL: 150 kg/330 lbs

**Direct Healthcare
Group**

Advancing Movement & Health®

**EN Visual inspection**

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EN - Please read the instructions before you start using this product

Purpose

This instruction applies to shower stool with round seat. The shower stool is mostly used in the bathroom. It can be used to sit on at the bathroom sink, but it can also be used as a shower seat. The stool is very suitable in places with little space.

Adjustment – see drawing, page 10

The shower stool is height adjustable in 10 steps. The adjustment takes place by detaching the snaplocks from the legs one by one and selecting the desired height. Be sure to adjust the four legs to the same height and to press the snaplocks firmly onto the tubes again.

Maintenance

To be wiped clean after use. Clean the stool with domestic cleaning materials without abrasive or corrosive properties. Solvents and chlorine should not be applied.

Caution

- Max load: 150 kg.
- Only use the shower stool as described in the paragraph "Purpose".
- Do not stand on the shower stool!
- Do not push or pull the shower stool over the floor while the user is sitting on it.
- Check the correct position of the snaplocks before using the stool.
- Check the condition and the correct position of the snaplocks frequently.
- For optimal safety it is recommended to check the ferrules frequently and to replace them whenever necessary.
- Ensure that all legs have the same height.

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
Linido

Direct Healthcare Group's Linido product line is a wide range of high quality safety products designed to provide ergonomic support, comfort and increased independence for the user in various bathing, showering and toileting situations.

Contact your local distributor if you have any questions about the product and its use.

See www.directhealthcaregroup.com for a complete list of distributors.

Always make sure that you have the right version of the manual.

The most recent editions of manuals are available for downloading from our website, www.directhealthcaregroup.com

Direct Healthcare Group

Advancing Movement & Health®

Direct Healthcare Group B.V.
Weteringweg 9,
2641 KM Pijnacker, The Netherlands
Tel: +31 (0)15 369 54 40
Fax: +31 (0)15 369 71 64
www.directhealthcaregroup.com
info.nl@directhealthcaregroup.com

Direct Healthcare Group Sverige AB
Torshamnsgatan 35,
SE-164 40 Kista, Sweden
Tel: +46 (0)8-557 62 200
Fax: +46 (0)8-557 62 299
www.directhealthcaregroup.com
info.se@directhealthcaregroup.com