

## Tri-pos<sup>3</sup> alternating air cushion



Tri-pos<sup>3</sup> alternating air cushion provides 24-hour Very High Risk care in combination with one of our air mattress systems.

It is supplied with or without the pump unit.

### Key features at a glance

- Adjustable to patient's weight for optimum pressure relief.
- It is powered by the same pump unit as all of our air mattresses. The added advantage of the interchangeable pumps is that the existing air mattress can be held inflated while the cushion is in use.
- In the case of breakdowns, the advantage of interchangeable pumps allows you to use the cushion with any of your existing stock of First Technicare pump units.
- Zip opening is set inboard and underneath of the cover, which ensures the cushion being used the correct way up. It minimises the risk of contamination: underneath zips avoid long upper side seams which risk ingress of fluids; fluids will not rise upwards.
- Power unit is both energy efficient and quiet in operation with lightweight and compact construction.
- Security: adjustable girth straps fasten cushion securely to a chair.
- Non-slip base stops the cushion from downward slide.
- two ways stretch covers are waterproof/ vapour-permeable for optimum patient comfort and care.
- **Technical Specifications**

### Pump

|                                    |  |
|------------------------------------|--|
| <b>Order Code</b>                  | TP3SPump IT                            |
| <b>Dimensions</b>                  | 305 x 145 x 100 mm                     |
| <b>Weight</b>                      | 1.8kg                                  |
| <b>Rated Voltage</b>               | 220/240V                               |
| <b>Cycle Time</b>                  | 15 minutes                             |
| <b>Alteration Cycle</b>            | ABAB                                   |
| <b>Electricity Safety Standard</b> | BS 5742, Part 1; IEC 601-1; EN 60601-1 |

### Cushion

|                                    |                                      |
|------------------------------------|--------------------------------------|
| <b>Order Code</b>                  | TPC (without pump), TPCP (with pump) |
| <b>Dimensions</b>                  | 530 x 508 x 76 mm                    |
| <b>Weight</b>                      | 1.5kg                                |
| <b>Maximum Patient Weight</b>      | 200kg (32 stones)                    |
| <b>British Standard Compliance</b> | BS6807; BS7175; BS3424/37            |