

# Direct Healthcare Group

Advancing Movement & Health®

## Rental & Service Solutions

### Equipment Guide

Bariatric Beds & Mattresses

Bariatric Seating & Cushions

Specialist Seating

Bariatric Wheelchairs

Walking Frames & Stand Aids

Hoisting & Handling

Scale Systems

**EMERGENCY BARIATRIC SALES & RENTAL**

**0800 879 9289**



**nightingale**

# Direct Healthcare Group

Advancing Movement & Health®

## WE ALL NEED TO MOVE TO LIVE.

HUMAN MOVEMENT IS FUNDAMENTAL TO OUR HEALTH AND QUALITY OF LIFE.

Direct Healthcare Group uniquely enables the healthcare challenges associated with reduced patient movement to be anticipated, managed, and in many cases, resolved to support care providers and patients.

**All leading to better health outcomes, accelerated recovery and easier daily living.**

### Direct Healthcare Services

#### Intelligent Pressure Care

A leader in the health application of pressure care technologies, Direct Healthcare Services designs and manufactures a comprehensive range of clinically proven, best value pressure care solutions including mattress systems, cushions and overlays.



#### Specialist Seating

Kirton specialist seating for patients with complex care needs are designed for comfort, support and stability. Incorporating clinically proven pressure care, Kirton seating solutions effectively promotes mobility and independence vital to maintaining quality of daily living.



#### Rental & Service Solutions

24/7/365 Rental service for standard, bariatric and specialist seating equipment.

Technical Service solutions to maintain equipment to required quality and regulatory standards.

# Why Choose Nightingale?

**Professional & Experienced**  
20 years Bariatric care experience  
across multiple sectors



**Fully Compliant to  
Regulatory Standards**  
Ensuring highest quality of care



**Best Value Solutions when  
you need them**  
Under 4 hours standard delivery



**Specialist Bariatric  
Education & Tools**  
Supporting the delivery of dignified, harm-free care



**Innovative Pressure Care  
Equipment**  
Award winning & clinically proven



**Complete Bariatric Solutions**  
Continuity of care across the  
full patient pathway



Beds



Mattresses



Seating

Specialist  
Seating

Cushions

Moving &  
Handling  
Equipment

Wheelchairs

Shower &  
BathroomWalking  
Frames/  
Stand Aids

Scales

## Equipment Selection Guide for Bariatric Patients

### MOBILE PATIENT



+

Olympia or Floor Bed  
with Mercury Advance Mattress

+

Static Chair or Pro Axis Riser Recliner  
Mercury Advance Air CushionCommode  
(choose suitable size)Walking Frame/  
Standing Aid

### SEMI MOBILE PATIENT

Gantry or Mobile Hoist  
with Lift Pants/Sling

Standing Aids

Return 7600 or  
Walking/Stand Aid

### NON MOBILE PATIENT



+



+

Olympia or Floor Bed with Mercury Advance /  
Bari Mattress & Dyna-Tek Heel Boot (Pair)Low Air Loss or  
Low Air Loss Turn AssistGantry/Mobile Hoist  
with SlingRepositioning  
SheetHigh Dependency  
ChairWheelchair Chair  
(choose option)

## BARIATRIC BED RANGE



### Olympia Bariatric Bed (500kg SWL)

#### Acute Bariatric Width-Adjustable Bed Frame

The Olympia Bariatric width-adjustable, profiling bed frame is fully compliant to the latest hospital bed standards. It can be used as a standard bed at 900mm and is expandable to 1100mm, 1200mm and 1300mm, with an extendable length to 2200mm.

The Olympia Bariatric has a user friendly, quiet operation to assist with the provision of dignified care.


**Order Code:** 0111

**Specifications:**

- 78st/500kg SWL (68st/435kg max patient weight)
- Width Adjustable from 900mm to 1300mm (with 1100 and 1200 options)
- Extendable length up to 2200mm

**Rental Daily Rate:** £44.00

**Profiling Positions**



### Raising Aid

#### Optimal support when getting in and out of the bed

Weight tested to be used as a handling aid.

Ergonomic solid laminate design.

**Order Code:** 1113 (Included with bed rental by request)



### Bariatric Floor Bed (300kg SWL)

The bariatric floor bed is a fully profiling bed frame offering the lowest bariatric frame height currently available of 100mm to greatly assist with patients that can be mobilised in and out of bed and for vulnerable patients that are prone to falls.

Clear access is also given under the bed frame for mobile hoists. The EN9000KS bed has a user friendly and quiet operation to assist with the provision of dignified care.

Comes with integral safety rails and removable head/footboards.

**Order Code:** 0112

**Specifications:**

- 47st/300kg SWL
- Frame width 1100mm

**Rental Daily Rate:** £44.00





## Pro Axis 2 (400kg & 500kg SWL)

### Expandable Bariatric Bed with Integrated Weigh Scale

The Pro Axis 2 is a 4-section (3 jointed, 1 fixed centrally), adjustable, profiling bed with integrated patient weighing system. The mattress platform can be adjusted to three widths, 1020mm, 1120mm and 1220mm, whilst also featuring a pull-out bed lengthener allowing for a total extension length of 1800mm.

As standard, the bed is equipped with a push-button handset for the patient and with an operator control handset for the exclusive use by specialised personnel with enables back rest adjustment, knee break adjustment, height adjustment, cardiac chair and massage position, plus Trendelenburg and reverse-Trendelenburg push-buttons.

The Precision Class III weighing system features four sensors, integrated into the structure of the bed, at the four corners to limit sensor deflection and stress, granting durability and better load distribution. The fixed support display, with rotating viewer, is battery powered allowing for a 4-hour operating time without a power connection.


**Order Code:** 0113

#### Specifications:


- 63st/400kg & 78st/500kg SWL
- Width Adjustable from 1020mm to 1220mm (with 1220 option)
- Extendable length up to 2180mm

**Rental Daily Rate:** £44.00

#### Profiling Positions



**Intelligent®**  
Pressure Care  
Management



## BARIATRIC MATTRESS RANGE

### Dyna-Form® Mercury Advance Bari WA (380kg Weight Limit)

A width-adjustable mattress replacement system that can deliver both pressure redistribution in static mode and pressure relief in dynamic mode, the Dyna-Form® Mercury Advance Bari WA is designed for use with patients at Very High Risk of developing pressure damage.

Offering high levels of patient comfort, this unique system has the facility to 'step up' to that of a dynamic mattress when clinically required. Similarly the mattress function can be 'stepped down' as the patient's condition improves. Can be used at 900mm, 1100mm, 1200mm and 1300mm widths.

**Mercury Advance Bari Mattress Order Code:** 0415

**Mercury Advance Bari Pump Order Code:** 1209

**Risk Category:** Very High Risk

**Dimensions:** 198cm (L) x 90–110/120–130cm (Width Adjustable) x 15cm (D)

**Rental Daily Rate (Mattress & Pump):** £22.00



### Dyna-Form® Low Air Loss (453kg Weight Limit)

The Dyna-Form® Low Air Loss system is an effective True Low Air Loss system designed for patients with sensitive skin and at Very High Risk of developing pressure ulcers.

Incorporating a high performance blower, the Dyna-Form® Low Air Loss system delivers true Low Air Loss therapy and excellent microclimate control with an evenly distributed flow of cooling air maintained across the mattress surface to assist in managing the heat and humidity of the skin.

**Low Air Loss Mattress Order Code:** 0403 **Low Air Loss Pump Order Code:** 1203

**Risk Category:** Very High Risk

**Dimensions:** 200cm (L) x 90/110/120/130cm (W) x 25cm (D)

**Rental Daily Rate (Mattress & Pump):** £22.00



### Dyna-Form® Power-Turn Elite Lateral Rotation Mattress (453kg Weight Limit)

Designed for challenging bariatric patients that have very little or no mobility, the Dyna-Form® Power-Turn Mattress System combines two clinically effective therapies: Lateral Rotation and True Low Air Loss in a portable system. A full, Left or Right Turn Mode is offered with either a 25° Small Turn or a 40° Big Turn.

**Order Code:** 0413

**Risk Category:** Very High Risk

**Dimensions:** 200cm (L) x 110cm\*/120cm\*\* (W) x 25cm (D)

**Rental Daily Rate (Mattress & Pump):** £27.50

**Important Note:** \*110cm (W) has a 363kg Weight Limit, \*\*120cm (W) has a 453kg Weight Limit

**Intelligent**  
Pressure Care  
Management



## BARIATRIC CUSHION RANGE

### Mercury Advance Bari Alternating Air Cushion

The Dyna-Form® Mercury Advance Bari Alternating Cushion is a dynamic replacement seating system, combined with the benefits of modern foam technology.

It is specifically designed for patients considered to be at Very High Risk of pressure ulcer development.

Offering high levels of patient comfort, this unique system has the facility to 'step up' to that of a dynamic cushion when clinically required. Similarly, the cushion function can be 'stepped down' as the patient's condition improves.

**25" seat width (318kg)**  
**Order Code:** 1101  
**Rental Daily Rate:** £5.50

**28" seat width (318kg)**  
**Order Code:** 1101  
**Rental Daily Rate:** £5.50

**30" seat width (318kg)**  
**Order Code:** 1101  
**Rental Daily Rate:** £5.50

**32" seat width (318kg)**  
**Order Code:** 1101  
**Rental Daily Rate:** £5.50

**34" seat width (318kg)**  
**Order Code:** 1101  
**Rental Daily Rate:** £5.50

## BARIATRIC SEATING RANGE



### Pro Axis Riser Recliner Chairs

**Hoist compatible – Option to hire integral Mercury Advance with riser recliner**

These heavy duty chairs are ergonomically designed and offer pressure relieving seats, footrests and waterfall backs with full electric rise and recline functionality. All chairs are also compatible with mobile hoists and stand aids. Adjustable seat depth up to 5" with dynamic back rest adjustment to accommodate different body shapes.

**25" seat width (318kg)**  
**Order Code:** 0504  
**Rental Daily Rate:** £27.50

**31" seat width (318kg)**  
**Order Code:** 0503  
**Rental Daily Rate:** £33.00

**28" seat width (318kg)**  
**Order Code:** 0504  
**Rental Daily Rate:** £27.50

**34" seat width (318kg)**  
**Order Code:** 0503  
**Rental Daily Rate:** £33.00





### Static Bedside Chair

Comes with drop arms and height adjustment

**28" seat width (318kg)**  
**Order Code:** 0502  
**Rental Daily Rate:** £16.50

**32" seat width (318kg)**  
**Order Code:** 0501  
**Rental Daily Rate:** £16.50



### Fero Obesitas Bedside Mobile Chair

A versatile bedside chair with electric backrest recline operation and drop arms. This chair can also be used as a patient transfer chair as it comes with larger castors, rear handle and a footplate.


**30" seat width (318kg)**  
**Order Code:** 0507  
**Rental Daily Rate:** £27.50



### Hi-Lo Tilt-In Space Bariatric Day Chair (300kg SWL)

Offering the capability and functionality of a bed with full electric operation and a high dependency castellated foam surface. This chair is designed to be used as the first step from bed to chair for high dependency patients and allows for easier bed to chair lateral transfers. Can be used to assist with sit to stand exercises during early stages of a patients rehabilitation. Large castors and a rear handle also allow for patient transportation.

**Order Code:** 0508  
**Rental Daily Rate:** £49.50



## COMMODOES/SHOWERING & WHEELCHAIR RANGE



The Freeway T100 Bariatric Shower Chair is suitable for a variety of applications for bariatric clients. The T100 benefits from increased diameter tubing, special heavy duty castors, heavy duty footrests with height adjustable heel supports which are removable and swing away.

Other commodes/showering accessories available.  
 Standard Commodes available from RDR: £10

**28" seat width (350kg)**  
**Order Code:** 0706  
**Rental Daily Rate:** £16.50

**30" seat width (350kg)**  
**Order Code:** 0706  
**Rental Daily Rate:** £16.50

**32" seat width (350kg)**  
**Order Code:** 0706  
**Rental Daily Rate:** £16.50





A range of folding and fixed width wheelchairs available to aid patient mobilisation.

**20" tilt in space wheelchair (125kg)**

**Order Code:** 0606

**Rental Daily Rate:** £22.00

**28" fixed width wheelchair (325kg)**

**Order Code:** 0606

**Rental Daily Rate:** £16.50

**30" folding wheelchair (318kg)**

**Order Code:** 0606

**Rental Daily Rate:** £16.50

**28" folding wheelchair (318kg)**

**Order Code:** 0606

**Rental Daily Rate:** £16.50



## ROLLATORS, WALKING FRAMES & STAND AIDS

### Stand Tall XXL Walking Aid (325kg SWL)

Height adjustable

**Order Code:** 0901

**Rental Daily Rate:** £16.50



### Height Adjustable Walking Aid (318kg SWL)

Height adjustable

**Order Code:** 0902

**Rental Daily Rate:** £5.50



### Return 7600 (205kg SWL)

Assisting with sit to stand  
and patient rehabilitation

**Order Code:** 0904

**Rental Daily Rate:** £11.00

## HOISTING & HANDLING



### Liko Ultra Twin Gantry Hoist (400kg SWL)

Includes re-usable sling

With a very high lifting capacity this gantry hoist has two lift motors that work together to meet very high demands. The Ultra Twin not only lifts and transfers, it also helps position the patient in sitting, semi-reclining or horizontal positions. Adjustable width span of 2m, 2.5m and 3m available.

**Order Code:** 0602 **Rental Daily Rate:** £49.50 **Additional Slings:** £3.85 per day



### Liko Viking XL Mobile Hoist (300kg SWL)

Includes re-usable sling

An increasing need for lifting heavy patients in different care situations led to the development of a mobile lift with high lifting capacity, made of aluminium to be as easy to manoeuvre as possible. This hoist comes with patient stability bars allowing for use as a rehabilitation frame with the Lift Pants as an accessory (please see accessories).

**Order Code:** 0601 **Rental Daily Rate:** £33.00 **Additional Slings:** £3.85 per day



### Titan 250 Bariatric Stand Aid (250kg SWL)

Includes re-usable sling

The Titan 250 is designed to carry out patient transfer from one seated position to another without compromising safety, for example in a toileting or general transfer situation. An ideal aid for patients who are cognitive and able to cooperate, and who are partially weight bearing but require assistance to be raised to an upright position, the Titan 250 facilitates an imitation of the body's natural rising movement, making the lifting action comfortable for the client.

**Order Code:** 0607 **Rental Daily Rate:** £33.00



### Liko Ultra Twist Plate (500kg SWL)

For use with Gantry Hoist only

Allowing for safe patient rehabilitation between the span of the gantry hoist system with 360 degree rotation.

**Order Code:** 1115 **Rental Daily Rate:** £16.50



### Liko Ultra Lift Pants (500kg SWL)

For use with Gantry or Viking XL Mobile Hoist

Lift Pants can be used with an overhead lift or with a mobile lift with sufficient lifting height allowing for safe patient rehabilitation.

Available in XXL or XXXL

**Order Code:** 0604 **Rental Daily Rate:** £3.85



### Liko Repositioning Sheet (500kg SWL)

For use with Gantry or Viking XL Mobile Hoist

The RepoSheet is placed under the patient and can function like an ordinary sheet during the entire care period. With minimal effort, caregivers are then able to easily reposition and turn the patient with the help of an overhead lift.

**Order Code:** 0603 **Rental Daily Rate:** £3.85



### Mobile Flex 2 Post Gantry Hoist (275kg SWL)

Includes re-usable sling

**Order Code:** 1115

**Rental Daily Rate:** £49.50



### Atlas 250 Mobile Lift (250kg SWL)

Includes re-usable sling

**Order Code:** 0601

**Rental Daily Rate:** £33.00

## CLIENT WEIGHING SERVICE RANGE

We can offer individual equipment assessments for beds, pressure care systems, seating, hoisting and patient/client mobilisation when accompanied by a Clinician.

In addition to this we can also offer a patient/client weigh service. (Please enquire for more information).



### Independent Bed Scale System (1000kg SWL)

#### Class III Compliant

A versatile low profile system allowing for patient weighs on virtually any bed frame. Comes with an integral trolley for ease of transportation.

Simple to use with BMI function.

**Order Code:** 0802

**Rental Daily Rate:** £16.50



### M200 Chair Scales

#### (300kg SWL)

#### Class III Compliant

Heavy duty wheeled chair scales with a high capacity wider seat, hinged armrests and foot plates. This chair allows for dignified patient weighs with a simple to use display unit which includes a BMI function.

**Order Code:** 0803

**Rental Daily Rate:** £11.00



### Ramped Foldaway Scale System (300kg SWL)

#### Class III Compliant

Suitable for wheelchair or standing weighs.

Transportation wheels and foldable for storage.

**Order Code:** 0804

**Rental Daily Rate:** £11.00

## TRAINING RANGE

### Bariatric Training Suit

**Order Code:** 1106

Provided to help support a dignified approach to bariatric care and assisting carers with training days.  
(Please enquire for more information).



## POSITIONING PRODUCT RANGE



The Dyna-Tek® Heel Boot and Dyna-Tek® Heel Pad are simple-to-use solutions that can help reduce the risk of pressure ulcers developing on the vulnerable heel area. Suitable for patients evaluated to be at High Risk of developing pressure damage, the Dyna-Tek® Heel Boot and Dyna-Tek® Heel Pad have been designed to offer patients comfortable and effective pressure redistribution.

Lightweight, soft and yet very durable, the Dyna-Tek® Heel Boot and Dyna-Tek® Heel Pad are also economic preventative solutions, especially when weighed against the significant costs that can be incurred in relation to pressure ulcer treatment.

The Dyna-Tek® Heel Boot and Dyna-Tek® Heel Pad are single patient use products that can remain in place for up to 28 days.

### 1. Dyna-Tek® Heel Boot (Large)

**Order Code:** DYNAHEE001 **Purchase Only:** £159.44 (1 box / 36 pairs)

### 2. Dyna-Tek® Heel Pad

**Order Code:** DYNAHEE002 **Purchase Only:** £77.32 (1 box / 18 units)

## ACUTE & COMMUNITY BED RANGE



### Olympia Hi-Lo Acute Care Bed (230kg SWL)

The Olympia Hospital XLow bariatric bed is a perfect symbiosis of safety and comfort. A low bed position of  $\pm 25$  cm reduces the risk of falling to an absolute minimum, whereas the ergonomic working height of  $\pm 80$  cm guarantees efficient care.

**Height adjustment:** 25 – 80cm

**Order Code:** 0112

**Rental Daily Rate:**

£13.75 (frame only)

£16.50 (with static mattress)

£22.00 (with dynamic mattress)



### Aerys Hi-Lo Community Bed (230kg SWL)

The Aerys Hi-Lo bariatric bed combines safety and comfort. A low bed position of  $\pm 24.5$  cm reduces the risk of falling to an absolute minimum, whereas the ergonomic working height of  $\pm 80$  cm guarantees efficient care. The Aerys bariatric bed has been especially designed to meet the needs and the requirements of Acute hospitals, Long-term care units and Nursing homes/Care homes and being easily transferable through standard doorways it can be installed where other beds could prove difficult.

**Height adjustment:** 25 – 80cm

**Order Code:** 0113

**Rental Daily Rate:**

£13.75 (frame only)

£16.050 (with static mattress)

£22.00 (with dynamic mattress)



### Acute Floor Bed (200kg SWL)

The EN9000KS bariatric floor bed is a fully profiling bed frame offering the lowest bariatric frame height currently available of 100mm to greatly assist with patients that can be mobilised in and out of bed and for vulnerable patients that are prone to falls.

Clear access is also given under the bed frame for mobile hoists. The EN9000KS bed has a user friendly and quiet operation to assist with the provision of dignified care. Comes with integral safety rails and removable head/footboards.

**Height adjustment:** 10 – 72.5cm

**Order Code:** 0114

**Rental Daily Rate:**

£22.00 (frame only)

£24.75 (with foam mattress)

£29.15 (with dynamic mattress)





## XXL Xpress Bariatric Care Bed (270kg SWL)

The XXL Xpress medical bariatric bed combines both the form of a modern design with the functionality needed to assist carers in giving superior levels of care to patients up to 40st/250kg in weight. It helps to reduce the risks associated with the manual handling of this patient group. The extra width provides a patient surface that is 4 feet in width or more (1200mm or 1400mm) in a bed that can be easily transported in and out of community environments.

**Order Code:** 0110

**Rental Daily Rate:** £22.00

Plus mattress option if required.

**Dimensions:** 200cm (L) extendable to 220cm (L) x 120cm (W)

**Intelligent<sup>®</sup>**  
Pressure Care  
Management



## STANDARD MATTRESS RANGE

### Dyna-Form<sup>®</sup> Mercury Advance (254kg Weight Limit)

The Dyna-Form Mercury Advance is a dynamic mattress replacement system combined with the benefits of modern foam technology. It is specifically designed for patients considered to be at 'Very High Risk' of pressure ulcer development. Offering high levels of patient comfort, this unique system has the facility to 'step up' to that of a dynamic mattress when clinically required. Similarly, the mattress function can be 'stepped down' as the patient's condition improves.

**Mercury Advance Mattress and Pump Order Code:** 0416

**Risk Category:** Very High Risk

**Dimensions:** 198cm (L) x 88cm (Width) x 15cm (D)

**Rental Daily Rate (Mattress & Pump):** £7.00

**Intelligent<sup>®</sup>**  
Pressure Care  
Management



## STANDARD CUSHION RANGE

### Dyna-Form<sup>®</sup> Mercury Advance (254kg Weight Limit)

The Dyna-Form<sup>®</sup> Mercury Advance Alternating Cushion is a dynamic replacement seating system, combined with the benefits of modern foam technology. It is specifically designed for patients considered to be at "Very High Risk" of pressure ulcer development.

Offering high levels of patient comfort, this unique system has the facility to "step up" to that of a dynamic cushion when clinically required. Similarly, the cushion function can be 'stepped down' as the patient's condition improves. A higher maximum weight capacity, up to 24 stone / 152kg, allows the product to meet the challenges of heavier patients.

**19" seat width (318kg)**

**Order Code:** 1101

**Rental Daily Rate:** £5.50

**22" seat width (318kg)**

**Order Code:** 1101

**Rental Daily Rate:** £5.50





## STANDARD SEATING RANGE

### Pro Axis Riser Recliner Chairs

**Hoist compatible – Option to hire integral Mercury Advance with riser recliner**

These heavy duty chairs are ergonomically designed and offer pressure relieving seats, footrests and waterfall backs with full electric rise and recline functionality. All chairs are also compatible with mobile hoists and stand aids. Adjustable seat depth up to 5" with dynamic back rest adjustment to accommodate different body shapes.

**19" seat width (120/140kg)**

**Order Code:** 0505

**Rental Daily Rate:** £22.00

**22" seat width (120/140kg)**

**Order Code:** 0505

**Rental Daily Rate:** £22.00

## SPECIALIST SEATING RANGE

All specialist seating manufactured by Kirton is now available through Nightingale, and comes with clinically proven pressure care cushions.



### Florien II (180kg SWL)

The Florien II is a dynamic and highly mobile seating system that enables greater independence by promoting correct posture and facilitating a supportive, comfortable seated position. Featuring manual tilt-in-space, leg rest elevation and removable side arms for easier hoist access and side transfers.

Available with lateral support for an extra £2 per day.

**Order Code:** 0707

**Rental Daily Rate:** £22.00



### G-2 (180kg SWL)

For the semi- and non-ambulant user with high posture and pressure care needs, the G-2 features electric tilt-in-space action and built-in, helical-motion lateral supports.

**400mm Order Code:** 0708

**450mm Order Code:** 0709

**500mm Order Code:** 0710

**550mm Order Code:** 0711

**Rental Daily Rate:** £22.00



### Duo Major (160kg SWL)

The Duo Major is designed for the non-ambulant user with profound posture and pressure care needs. Electric tilt-in-space and dynamic features provide extensive positioning and postural management for individuals who are at risk of developing an unconventional body shape, unable to maintain a midline position and require additional support.

Available with lateral support for an extra £2 per day.

**350mm Order Code:** 0712

**400mm Order Code:** 0713

**450mm Order Code:** 0714

**500mm Order Code:** 0715


**Rental Daily Rate:** £22.00

# Equipment Selection Guide


Our professional teams are always on hand to help you choose the best solution for your needs. Here are some patient details you might want to consider when selecting bariatric equipment:

## Body Shape

Apple



Pear Abducted



Pear Adducted



Bulbous Gluteal



Age Group	Child	Young Adult	Adult	Mature Adult	Elderly
Height	to 5ft	5ft to 5ft 6"	5ft 6" to 6ft	6ft to 6 ft 6"	over 6ft 6"
Weight	to 20st	20–30st	30–40st	40–50st	50–65st

Mobility Assessment	Pressure Ulcer Risk Assessment
No Limitation (Can stand and walk)	Low
Slightly Limited (Walks with support)	Low to Medium
Limited (Can bear weight)	Medium
Very Limited (Chair bound)	Medium to High
Completely Immobile (Bed bound)	High
Not Conscious	High to Very High Risk

## Simple Patient Assessments make equipment choice easier and more cost effective

Order Maker: (Name and Title)	Tel:
Hospital/Ward:	Email:
Patient Name/ID:	Weight:                      Height:
Body Shape:	Mobility:
PAC Needs:	Date and Time:
Ward Contact:	Tel:
Trust:	ORDER NO:

## Always order equipment in line with local protocol and clinical assessment

Notes:





# Direct Healthcare Group

Advancing Movement & Health®

Intelligent Pressure Care

Specialist Seating

Rental & Service Solutions



## Direct Healthcare Group

Withey Court, Western Industrial Estate  
Caerphilly, United Kingdom  
CF83 1BF

**T:** +44 (0) 800 043 0881

**F:** +44 (0) 845 459 9832

**E:** [info@directhealthcaregroup.com](mailto:info@directhealthcaregroup.com)



LIT-00041 Issue 1  
Date: January 2018

[DIRECTHEALTHCAREGROUP.COM](http://DIRECTHEALTHCAREGROUP.COM)

