



Chatsworth Advance®

The Chatsworth Advance® riser recliner is the ultimate bespoke seating solution. It comes with a range of options to suit your requirements and budget.

Choose from a single, dual or dual tilt-in-space action, each providing a smooth, effortless rise and recline movement, operated by an easy to use handset. The Chatsworth Advance® has the ability to adjust its lift and tilt for the best rising position to suit each user. It is supplied as standard in maximum tilt, lowest lift position but may easily be adjusted to achieve higher flatter lift or very high virtually flat lift upon request. This can be helpful with side transfers into bed or another chair.

The Chatsworth Advance® can be tailored to each individual's specification in terms of seat width, height, depth and armrest height to ensure comfort and support when seated and a choice of traditional scrolled or contemporary armrests to match the chair to your own style.

Technical Specifications

SWL	158kg/25 stone
Seat Height	440mm, 470mm, 495mm
Seat Width	470mm, 520mm, 520mm
Seat Depth (Single Motor)	445mm, 470mm, 495mm
Seat Depth (Dual T-i-S Motor)	495mm, 470mm, 495mm
Arm Height	190mm (single motor), 200mm (dual motor, dual T-i-S)
Back Height	650mm, 700mm, 700mm

Features

Benefits

Choice of motor

Gives the user a choice of functionality

Four backrest options

Provides optimal comfort and support

Bespoke size options

Creates a tailored chair for the users needs

Choice of pressure care options

Ensures the users pressure risks are met