



# Dyna-Form® Community

The **Dyna-Form® Community** mattress can be used with patients considered to be up to that of “**High Risk**” of pressure ulcer development. Using modern foam designs, castellated on one side to aid pressure ulcer management, whilst providing excellent patient comfort; this mattress can be used without a turning programme, aiding moving & handling procedures. The **Dyna-Form® Community** mattress can work comfortably with any electric profiling bed frame. Its durable design uses only the very best CME foam & fabrics to provide a modern mattress complete with infection control enhancing cover.

Also Available in Pro-Formance® Strikethrough Resistant Cover.

## Features

## Benefits

1a 1b	Castellated foam (with tear duct cut)	1a smaller castellated cells for heels and 1b larger cells sacrally
2	Single sided cut surface	Works without a turning programme
3	Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear & friction forces, whilst assisting infection control policies
4	Fully sewn or high frequency welded cover. (NB: welded cover option carries a small price premium)	Prevents fluid ingress & aids decontamination processes (can be laundered up to 95°C)
5	Fully replaceable parts	Ensures great value for money

## Technical Specifications

<b>Product Code</b>	MAT2410001
<b>Risk Category</b>	High Risk
<b>Weight Limit</b>	22 Stone / 140kg
<b>Warranty</b>	1 Year
<b>Dimensions</b>	(Length) 198cm x (Width) 88cm x (Height) 15cm – However Variable Sizes Available Upon Request
<b>Foam Density</b>	33kg/m <sup>3</sup> Nominal Hardness – 130 Newtons
<b>Fatigue Class</b>	Very Severe
<b>Fire Retardancy</b>	BS7177; 2008 (Crib 5 Composite)
<b>Weight</b>	10kg
<b>Colour</b>	Blue
<b>Hypochlorite Solution compatible</b>	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)