


FlexibleSling

SystemRoMedic™

Trial Fitting Protocol

Max: 300kg



Patient's name: _____

Social insurance number: _____

FlexibleSling is pliable and flexible around the body and provides a feeling of security. It is available in three models; one with divided leg supports, one with undivided leg supports and one with undivided leg supports and commode opening. FlexibleSling is designed to give a lot of support and to adapt to the shape of the patient's body. After use, it can be left in the wheelchair, moulded seat or similar.



Functional inspection

Visual and mechanical inspection

Check the condition and function of the sling regularly. Always inspect the product after laundering. Check to ensure that seams and material are free from damage. Check to ensure that fabric is not worn or faded. Apply load to the device and check to ensure that clasps, handles, etc. withstand heavy load. If there are signs of wear, the product must be discarded.

Which size is used?

XX-Small X-Small Small Medium Large X-Large XX-Large

Which material is used?

Polyester net

Accessories

SoftLegSupport AdjustmentLoop StrapsPadding

Which slingbar is used?

Leg support placement (FlexibleSling with divided leg supports):



Crossed



Overlapped



Each leg separate

Instructions for use

Can the sling be left under the patient in the wheelchair? Yes No

Can the sling be used for transfers to/from the toilet? Yes No

Can the sling be used in showering situations? Yes No

Loop alternatives/AdjustmentLoop: _____

Have staff that use the assistive device:

Read the instructions? Yes No

Received practical instruction in patient transfers? Yes No

Instructions to personnel:

e.g.; Raise the bed side rail on the side towards which the patient is to be turned. Placement of head/arms/hands?
Placement of legs/feet? Use bed functions.

Number of caregivers: _____

Instructions to patient:

e.g.; In what way can the patient participate actively in the transfer situation? Which instructions are used?

Other:

Trial-fitted by:

Date: _____

Phone: _____

Phone: _____

In the events of changes or questions, contact the trial-fitter.

Manuals and instructions for functional testing are available on DHG's website:

DHG-HEALTHCARE.COM

SystemRoMedic™ is a complete solution for easy transfers.
The concept encompasses four product areas that meet
most needs in different transfer situations.

Transfer, between two locations. **Positioning**, within the same location.
Support, for mobility. **Lifting**, both manual and mechanical.

DHG
Moving Health Forward

Withey Court, Western Industrial Estate
Caerphilly, United Kingdom CF83 1BF
T: +44 (0) 800 043 0881
F: +44 (0) 845 459 9832
E: info@directhealthcaregroup.com
DHG-HEALTHCARE.COM