




Mover Aqua

Mover is a transfer platform developed to provide smooth, convenient transfer transfer indoor/flat surfaces. It's easy to use and helps transfer patients safely and conveniently. Mover has a narrow base frame and double castors for easy manoeuvrability even in confined spaces. Mover has a stable platform and comfortable seat cushions, and thanks to constant eye contact with the carer, the user always feels comfortable and confident. The new Mover Aqua is specially designed to be used also in showers. PU cushions together with stainless solutions makes this possible. Hygiene management is improved as well.



Technical specifications

Part No.	Product	Max user weight	Width - external	Width - step	Height	Length	Seat height	Weight
53-301	Mover Aqua	150kg	63.5cm	38cm	107cm	94cm	83cm	28kg