

Mover Aqua

Mover Aqua is a transfer platform developed to provide smooth, convenient transfer on indoor/flat surfaces and has been specifically designed to be used also in showers, reducing the amount of equipment or transfers required.

SWL
150kg/
23 stone

Key features

Promotes **EARLIER MOBILISATION** by encouraging patients to be active during sitting to standing transfers

EASY TO USE and helps transfer patients safely and conveniently

The **ERGONOMICALLY DESIGNED HANDLEBARS** make steering comfortable and convenient for the carer

3 CROSSBARS at different heights allow the user to find their own individual grip, for safer standing

The **SHIN SUPPORTS** have large contact surfaces, which makes standing up comfortable for users no matter how long their lower legs. The shin supports are also padded to prevent pressure on the lower legs


The **NON-SLIP PLATFORM** surface helps users stand steadily without shoes

Stable platform and **COMFORTABLE SEAT CUSHIONS** – constant eye contact with the carer ensures the user always feels comfortable and confident.

The **PIVOTING SEATS** can be moved out to enable the patient to stand, and provides an angled seated support during the transfer

NARROW BASE FRAME and double castors for easy manoeuvrability, even in confined spaces

The **LOW STEP** means users need not lift their feet too high when getting ready to stand up



Technical specifications

Max. User Weight	150kg
Width – External	63.5cm
Width – Step	38cm
Height	107cm
Length	96cm
Seat Height	80cm
Weight	28kg