

LECKEY[®] Pal

Classroom Seat



Posture + Attention = Learning

Designed by the Leckey team of experienced therapists and engineers, Pal is an attractive, robust and easy to use classroom and nursery seating system for kids of 1-8 years of age, with mild to moderate postural needs. It is designed to give improved stability, which helps reduce fatigue and allows the child longer periods of concentration and fine motor activity. Available in 3 sizes and 4 colours it will help kids up to 8 years of age to participate, learn and give their full attention in the classroom.






Created for self assessment
 The simplicity of Pal enables the therapist to assess the child's needs, instantly order their Pal Classroom Seat and quickly and easily set up the product when it arrives. Pal really is the simple, fast way to help kids join in with learning.

Ordering is as easy as 1 2 3

1. Measure
 Pal is a self assessment product where you assess the size required by 3 simple measurements:

1. Seat depth
2. Seat width
3. Seat height



Pal is available in 3 sizes and comes in 4 fab colours: green, orange, blue and pink.

2. Order
 Complete the Pal order form indicating:

1. Chair size and colour.
2. Any leg accessory: Castor, Rocker, Stabiliser or High Seat.
3. Any other accessory items: Footrest, Tray, Table, Pommel or Thick Side Pads.

3. Delivery
 Our fast track service will deliver your chair within 14 days. The chair will be delivered assembled and ready for use.



Pal Seat shown with footrest accessory



Pal Seat shown with castor leg accessory



Pal Seat with rocker leg accessory



Pal Seat with stabiliser leg accessory



Pal Seat with high leg accessory (Size 1 only)

Pal Sizing Chart

	Age Range	Stature Range	User Weight	Seat Depth	Seat Width	Seat Height	Backrest Height	Seat Weight	Table Height
Size 1	1 - 3 years	700 - 950mm	Max - 30kg	180 - 240mm	150 - 250mm	160 - 240mm	265mm	4.3kgs	400mm
Size 2	2 - 5 years	900 - 1050mm	Max - 35kg	220 - 290mm	170 - 270mm	185 - 285mm	315mm	6kgs	460mm
Size 3	4 - 8 years	1000 - 1300mm	Max - 40kg	250 - 330mm	230 - 330mm	250 - 350mm	380mm	8.2kgs	530mm


Backrest Angle 95 - 105°

The Pal Seat comes as standard with depth adjustable seat base, angle adjustable back rest, 2 point pelvic harness, height adjustable legs and wipe clean back rest and seat base cushions.

Pal Product Codes

X = Size 1, 2 or 3


AA = Colour (01 Green, 02 Orange, 03 Blue, 04 Pink)




Seat (including Lap Strap)
135-X600-AA



Ski with Castors (Pair)
135-X610



Rocker (Pair)
135-X604



Stabilisers (Pair)
135-X604



High Seat Frame (Size 1 only)
135-1630



Footrest
135-X603



25mm Side Pad (Pair)
135-X635-AA
50mm Side Pad (Pair)
135-X654-AA




Activity Tray
135-X620
Tray Pad (Not shown)
135-X656



Pommel
Size 1 - 135-1602-AA
Size 2 & 3 - 135-3602-AA




Ramped Seat Cushion
135-X651-AA



Tables are also available in the four colours: green, orange, blue and pink.

4 Person Table
1300x500mm
135-X701-AA



2 Person Table
650x500mm
135-X700-AA

LECKEY[®]

Leckey Healthcare
Kilwee Business Park
Dunmurry
BT17 0HD
Northern Ireland
United Kingdom

T: 028 9060 2277
F: 028 9060 0795
E: info@leckey.com
W: www.leckey.com

LS210-01

