

QuickLink®

SystemRoMedic™



Easy-to-use and versatile

QuickLink® provides an easy-to-use and versatile solution for the rapid changing of sling bars and associated accessories on mobile and ceiling lifts. The spring-loaded pull tab opens the system and securely closes it once the sling bar or accessory is in place, to ensure a secure connection.

Audible and tactile feedback

An audible click and tactile spring closure feedback provides the user with a clear indication that the hoist sling bar or accessory is safely attached.

Robust and assured design

A robust design and material provides the user with confidence that the hoisted user is safe and secure. A 300kg / 660lbs maximum user weight provides a large lifting capacity for use in a wide range of patient demographics. Testing to ISO 10535:2006 provides assurance that QuickLink® is a proven design.

Technical Information

Article name	QuickLink®
Article number	70200101
Maximum user weight	300kg / 660lbs
Material	Forged, high quality steel
Weight	Approx. 420g
Testing standard	ISO 10535:2006
Expected lifetime	Minimum 20,000 lifts (approx. 5 lifts/day for 10 years)
Dimensions	60 x 82 x 30mm
Height affecting lift height	50mm
Compatibility, lifts	Vega, Carina, Eva, EvaDrive, RiseBasic, RiseAtlas
Compatibility, slingbars	SlingBar S, M, L 4-point SlingBar StretcherBar
Compatibility, accessories	Charder scale