



Backrest for shower seat

Code **G12JDS06** Series **HUG**

Net weight kg 1,6 Dimensions [WxDxH] mm 420x40x270

In use weight capacity Vertical kg

Use and installation info

Use a soft, moist cloth soaked with liquid soap and/or other liquid sanitizing substances. **DO NOT USE** abrasive liquid detergents, chlorine bleach or strong acids that does not improve the aseptic aspect of the equipment and could lead to opacity and deterioration of the product.

Fasteners

Ponte Giulio suggest to use fasteners suitable to the wall characteristics on which the product will be installed.

Packaging information

Dimensions [WxDxH] mm 395x130x450 Weight kg 2,05

Available colours

WE Signal White, **DD** Graphite, **PB** Coral