





Seat, with backrest and armrest

Code G12JDS04 Series HUG

Net weight kg 7,2 Dimensions [WxDxH] mm 536x527x750

In use weight capacity Vertical kg

## Use and installation info

Use a soft, moist cloth soaked with liquid soap and/or other liquid sanitizing substances. DO NOT USE abrasive liquid detergents, chlorine bleach or strong acids that does not improve the aseptic aspect of the equipment and could lead to opacity and deterioration of the product.

Fasteners

## **Packaging information**

Dimensions [WxDxH] mm 435x145x485 Weight kg 8,1

## **Available colours**

Primary: WE Signal White, DD Graphite, PB Coral

Secondary: G9 Lime, R2 Copper, P5 Aubergine, WE Signal White, N2

**Bright** 



