



Wendover Advance®

The Wendover Advance® riser recliner is a contemporary chair designed with simplicity in mind. Available in small, medium and large sizes, this elegant chair accommodates the majority of dimensions to achieve a comfortable seated position that promotes good posture.

Available with Soft/Cushioned Back, Waterfall Back, Lateral Support Back or T-Back to tailor the comfort level to suit individual preferences. The Wendover Advance® can be underpinned by a choice of action depending on your functionality needs; single motor for those who prefer to maintain a seat to back angle while reclining or dual tilt-in-space for the added function of an independent legrest and backrest movement.

Features

Benefits

Choice of motor and chair size

Creates a tailored chair for the users size and needs

Four backrest options

Provides optimal comfort and support

Contemporary armrest

Fits into most domestic environments

Choice of pressure care options

Ensures the users pressure risks are met

Technical Specifications

SWL	158kg/25 stone
Seat Height	440mm, 470mm, 495mm
Seat Width	470mm, 520mm, 520mm
Seat Depth (Single Motor)	445mm, 470mm, 495mm
Seat Depth (Dual T-i-S Motor)	495mm, 470mm, 495mm
Arm Height	190mm
Back Height	650mm, 700mm, 700mm