

Description	DHS Order Code	NHSSC Order Code	PIP Order Code
Dyna-Tek Heel Boot	HEEL/LIFT/PAIR Order Quantity: 1 Box (36 Pairs)	FER13763 Order Quantity: 1 Box (36 Pairs)	399-2369 Order Quantity: 1 Pair
Dyna-Tek Heel Pad	HEEL/PAD/PAIR Order Quantity: 1 Box (18 Pairs)	FER13762 Order Quantity: 1 Box (18 Pairs)	402-9096 Order Quantity: 1 Pair



Dyna-Tek™ Positioning Products

The **Dyna-Tek Heel Boot** and **Foot Protector** are simple-to-use solutions that can help reduce the risk of pressure ulcers developing on the vulnerable heel area. Suitable for patients evaluated to be at **'High Risk'** of developing pressure damage, the Dyna-Tek Heel Boot and Foot Protector has been designed to offer patients comfortable and effective pressure redistribution.

Lightweight, soft and yet very durable, the Dyna-Tek Heel Boot and Foot Protector is also an economical preventative solution, especially when weighed against the significant costs that can be incurred in relation to pressure ulcer treatment. The Dyna-Tek Heel Boot and Foot Protector is a single patient use product that can remain in place for up to 28 days.

Features

Durable castellated foam design

Easy to place and adjust
Velcro fastening

Lightweight soft foam cradles
and protects the foot and ankle

Crib 5 Fire Tested

Single patient use product
that can remain in place for
up to 28 days

Benefits

Effectively redistributes pressure by
up to 50%* and reduces the risk of
pressure damage

Quick and simple to use or place
and remove

Full positional support and maximum
patient comfort

Complies fully with the latest fire
testing regulations

A reliable and economical solution
that can be utilised as part of a
pressure ulcer prevention strategy

Technical Specifications

Product Code 1. **HEEL/LIFT/PAIR (36 pairs per carton)**
2. **HEEL/PAD/PAIR (18 units per carton)**

Risk Category High Risk

Fire Retardancy Crib 5

*Direct Healthcare Services, 2014, Pressure Mapping held on file