



Design – Patent Pending

Dyna-Tek™
Posture Visco

The **Dyna-Tek Posture Visco** Cushion is designed to create better postural management, along with great pressure re-distribution. This product is specifically made for users considered to be **“High Risk”** of pressure ulcer development. Equally at home in a wheelchair, ward or homecare environment, **Dyna-Tek Posture Visco** provides the ultimate combination of pressure care and patient comfort.

Features

- ① Uniquely sculpted “ischial tuberosity well” design
- ② Visco elastic foam surface
- ③ Finest quality CME foam
- ④ Multi-stretch, vapour permeable cover (impregnated with an anti fungal agent)
- ⑤ Non-slip base
- ⑥ Fully replaceable parts

Benefits

- Creates exceptional postural management
- Provides excellent pressure re-distribution & important patient comfort
- Prevents patient bottoming out and aids patient transfer
- Reduces shear & friction forces, whilst assisting infection control policies (can be laundered up to 95 °C)
- Aids moving & handling procedures, providing ultimate patient safety
- Ensures great value for money

Technical Specifications

Product Code	CUSH/PV/46/46/7.5
Risk Category	High Risk
Weight Limit	24 Stone / 152kg
Warranty	2 Years
Dimensions	(Length) 46cm x (Width) 46cm x (Height) 7.5cm (18” x 18” x 3”) Variable Sizes Available Upon Request
Foam Density	CME Base—38—40kg/m3 Nominal Hardness—200 Newtons
Fatigue Class	Severe
Fire Retardancy	BS7176; 2007
Weight	1.3kg
Colour	Blue
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)