

# Coral

Shower stool with  
back and arm supports

*Stable support*

*Height-adjustable*

*Compact and light-weight*

*Ergonomic design*

3002010130  
Coral, shower stool with back and arm supports



## **A compact and height-adjustable shower stool for independent comfort and convenience**

Coral shower stool with back and arm supports is an assistive device allowing the user to sit down safely and conveniently during showering or at the wash-basin. The shower stool provides a stable and height-adjustable support for active users with weak legs, low stamina and/or poor balance and who require a little more support during sit-to-stand. For maximum stability, the legs of the stool are fitted with non-slip end caps with drainage holes. The seat height is easily adjusted to ensure the perfect fit for every user and situation.



# Coral

## Molded shower seat provides a good sitting position and a comfortable seating angle

The rust-free frame is fitted with a practical and comfortable molded shower seat made from durable plastic with a smooth surface. The seat provides a good sitting position and the comfortable seating angle of one degree creates a safe situation and a pleasant experience during showering as the user will not slide off the seat.

## Back and arm supports provide extra support and security

For users with reduced strength and balance, sit-to-stand and stand-to-sit procedures can be difficult. This model of Coral features a firm, comfortable back support and fixed arm supports that provide the extra support needed during these movements. The back and arm supports also offer the user an extra sense of security while seated without hindering a more active and independent user.



## Light-weight and compact, easy to fit

The shower stool is light-weight and compact and fits easily into any wet room, also when space is limited. The design, with an opening in the back support, and the shape of the arm supports makes it easy for users with reduced strength or hand function to grip and to lift the stool.

## Also available with toilet seat

This model of Coral is also available with a toilet seat (item no. 3002010030). Coral stools are also available without back and arm supports.



For more information, please visit our website, [www.handicare.com](http://www.handicare.com), or contact your local Handicare partner.

## PRODUCT SPECIFICATIONS

Art.no	Description
3002010130	Coral, shower stool with back and arm supports

## FEATURES AND BENEFITS

- Stable and secure, ergonomic design
- Safe and very easy to use
- Compact, hygienic and attractive design
- Fits into any wet room, also when space is limited
- Very easy to assemble – no tools required
- Stainless steel and plastic, suitable for wet environments
- Easy maintenance
- Molded shower seat provides a good sitting position and a comfortable seating angle
- Easily adjustable seat height; 420-570 mm, adjustable in 7 steps of 25 mm
- Firm yet flexible and comfortable back support
- Fixed arm supports for support during transfers and for security while sitting
- Back and arm supports are designed not to hinder a more active and independent user

## TECHNICAL INFORMATION

Material, frame and legs: \_\_\_\_\_ Stainless steel  
Material, seat, back and arm supports: \_\_\_\_\_ GF 30 % PP plastic  
Seat height: \_\_\_\_\_ 420 - 570 mm / 16.5" - 22.4", adjustable  
Seat measurements, width x depth: \_\_\_\_\_ 490 mm x 410 mm / 19.3" x 16.1"  
Distance between arm supports: \_\_\_\_\_ 450 mm / 17.7"  
Height, back support: \_\_\_\_\_ 330 mm / 13"  
Width, back support: \_\_\_\_\_ 349 mm / 13.7"  
Height, arm supports (measured from seat): \_\_\_\_\_ 210 mm / 8.3"  
Length, arm supports: \_\_\_\_\_ 249 mm / 9.8"  
Total height: \_\_\_\_\_ 750 - 900 mm / 29.5" - 35.4", adjustable  
Outer measurements, width x depth: \_\_\_\_\_ 550 mm x 525 mm / 21.7" x 20.7"  
Weight: \_\_\_\_\_ 7 kg / 15 lbs  
Max. load on arm support: \_\_\_\_\_ 65 kg / 143.3 lbs per arm support  
Max. user weight: \_\_\_\_\_ 150 kg / 330 lbs

## ACCESSORIES

Item No.	Product
3002920501	Transportation seat, PUR, with cut-outs for fixed arm supports
3002920200	Calf support, textile, with Velcro



**handicare**

Maskinvägen 17,  
972 54 Luleå, Sweden  
Tel: +46 (0)8-557 62 200  
Fax: +46 (0)8-557 62 299  
[www.handicare.com](http://www.handicare.com)